



2023 Liberty Mutual Insurance Industrial Hygiene Lab sampling guide

AIHA LAP, LLC Laboratory ID # LAP-100045

Contact us for assistance:

lmihlaboratory@libertymutual.com

or

1-800-230-6263



Liberty Mutual Insurance Industrial Hygiene Laboratory (LMIHL) provides analytical services related to occupational disease exposure in their workplace.

Comprehensive service includes:

- Method development or validation
- Technical consultation
- Training
- Sample Analyses:
 - Aldehydes
 - Organic and Inorganic acids
 - Amines
 - Pharmaceuticals
 - Asbestos fiber counting and bulk identification
 - Silica
 - Medical gas analysis
 - Common and specialty metals
 - Grade D breathing air analysis
 - Common and specialty organic solvents
 - Scans – metals, solvent, acids, aldehydes, isocyanates, etc.
 - Gravimetrics for total, respirable or inhalable dust

Liberty Mutual SafetyNet™ is your source for chain of custody and media order forms as well as the Fee Schedule. You may also contact the lab.

Telephone assistance: 1-800-230-6263

Telephone assistance is available to all our customers from 7:30AM – 6PM, Eastern Time. We encourage our customers to call us with their technical questions relating to proper sampling and media selection. Listed below is our support team.

- Laura Melton, MS, Laboratory Manager, 508-544-5348
- Eva Longo, MS, Prod/QC Chemist Consultant, 508-544-5370
- Neela Joshi, MS Prod/QC Chemist Consultant, 508-544-5360

Analysis turnaround time

The Laboratory's standard turnaround time is 5 business days or better from the receipt of samples. Samples received after 3PM are considered next day's samples. Results are provided by email only. Rush analysis for routine samples can be arranged with prior notification.

The following surcharges apply for RUSH requests:

One business day	200% surcharge
Two business days	150% surcharge
Three business days	75% surcharge

Note: Certain Specialty Analysis cannot be done on a rush basis. Contact the laboratory at least 24 hours before shipping your samples for RUSH analysis.

Discounts for multiple analytes on the same media

The price listed for each individual analyte is the price for a single analyte per sample. Analysis for multiple compatible analytes collected on the same media is available at a discounted fee. The higher priced analyte is billed at the published fee with additional analytes discounted. Refer to the fee schedule.

Payment terms

The total cost of the services provided by the LMIH Laboratory will be based on the quoted rate. Prices are subject to change without notice. Payment term is net 30 days if you are paying by check.

Blanks

Blanks are required by all analytical methods and are used to measure contamination and media interferences. The recommended number of field blanks per sample set is 10% of the total number of samples, or a minimum of 1 blank per sample set. Blanks are analyzed following the same procedure as samples and are charged at the same rate.

Media/Shipping

Common media such as charcoal tubes, pre-weighed filters, etc. are provided at no extra cost as long as the samples and unused media are returned to the laboratory for analysis within 30 days. There is a media charge for specialty filters and sampling devices such as PPI and IOM. Unreturned media after 30 days will be charged at cost. Specialty media that cannot be reused, such as isocyanate filters, aldehyde badges and tubes, ozone filters, etc. will be invoiced at cost. Media are shipped UPS ground at no charge to the client. Rush and international shipments will be billed to the client at cost. Rush media orders for same day shipping must be submitted by 2 p.m. Eastern Time, Monday through Friday. Media requiring cold shipping cannot be shipped on a Friday.

- Request by email using the media order form.
- Return shipping is the responsibility of the customer.

Sample collection and sample submission

Sampling collection supplies

We provide free collection supplies (templates, gloves, plastic bags, coolers, etc.) and ship them to you via UPS ground. Rush and international shipments will be billed to the client at cost. It is the client's responsibility to ensure proper sampling, handling, packing (returning the cooler with frozen packs) and return shipping to the laboratory.

Chain-of-custody forms will be provided with every media and equipment order; please complete the form and submit with your samples. The chain-of-custody forms are also available on Liberty Mutual SafetyNet™.

Sample minimum

When preparing samples for analysis, please note that some analyses require a three (3) sample minimum fee to cover the cost of method set-up for non-routine analysis.

Special handling and shipping

Many NIOSH and OSHA methods indicate that solvent samples are stable at room temperature. LMIH Laboratory recommends shipping samples cold overnight especially during the summer months. Please refer to the list of Analytes for additional instructions. These guidelines must be followed in order to ensure the integrity and validity of the samples.

Bulk sample submissions

Certain analysis requires the submission of bulks. Bulks should always be packaged and shipped separately from the samples. Safety Data Sheets (SDS) must accompany the bulk samples.

Analytical results

LMIHL reserves the right to determine appropriate format in which the analytical results are reported. All results are provided for the exclusive use of our client. LMIHL accepts no responsibility or liability for the client's use of the analytical results.

LMIHL requires precise and complete instruction before it releases any reports. Any later request will require written permission from submitter.

Free pump loan program

We loan pre-calibrated sampling pumps and accessories to customers for a period of two (2) weeks at no cost if the samples are returned to the lab for analysis. All equipment is loaned on a first-come, first-serve basis. As a result, there may be instances when the equipment is not available at the time of your request. We suggest submitting your request for equipment at least one week prior to your scheduled sampling activity.

At the end of the two week period, all equipment must be returned and collected samples sent for analysis to the LMIH Lab. Equipment not returned after the two week period will incur an equipment rental charge. Contact your local Loss Control Representative or the LMIH lab for information/questions (lmihlaboratory@libertymutual.com or 800-230-6263) regarding your sampling needs.

Scans

The scans performed by the laboratory are listed below. We can also customize scans.

Aldehyde scan (Glutaraldehyde is collected and analyzed separately)			
Acetaldehyde	Benzaldehyde	Formaldehyde	Valeraldehyde
Acrolein	n-Butyraldehyde	Propionaldehyde	
Inorganic acid scan			
Bromide	Fluoride	Phosphate	
Chloride	Nitrate	Sulfate	
Isocyanate scan			
Hexamethylene diisocyanate (HDI)	Isophorone diisocyanate (IPDI)		
4,4-Methylene bisphenyl isocyanate (MDI)	2,4-Toluene diisocyanate (2,4-TDI)		
2,6-Toluene diisocyanate (2,6-TDI)			
Metals – Liberty's 20-Metal scan			
Aluminum	Cadmium	Lead	Tin
Antimony	Chromium	Magnesium	Titanium
Arsenic	Cobalt	Manganese	Thallium
Beryllium	Copper	Nickel	Vanadium pentoxide as V
Calcium	Iron	Selenium	Zinc
Metals – 13-Metal scan (welding fume scan)			
Antimony	Cobalt	Manganese	Vanadium pentoxide as V
Beryllium	Copper	Molybdenum	Zinc
Cadmium	Iron	Nickel	
Chromium	Lead		
Organic acid scan			
Acetic acid	Butyric acid	Formic acid	Propionic acid
Organic solvents (GC/MS) scan			
(Qualitative analysis. For quantitative analysis please call lab.)			
Acetone	Dimethyl formamide (DMF)	Methyl amyl ketone (MAK)	
Acetonitrile	1,4-Dioxane	Methyl cyclopentane	
Acrylonitrile	Epichlorohydrin	Methyl ethyl ketone (MEK)	
n-Amyl acetate	Ethanol	Methyl isoamyl ketone (MIAK)	

Benzene	Ethyl acetate	Methyl isobutyl ketone (MIBK)
Sec-Butanol	Propylene glycol methyl ether acetate (PGMEA)	Methyl propyl ketone (MPK)
n-Butyl acetate		Methyl t-butyl ether (MTBE)
s-Butyl acetate	Propylene glycol monomethyl ether (PGME)	a-Methyl styrene
n-Butyl acrylate		Methylene chloride
n-Butyl alcohol	Styrene	Methyl methacrylate
Butyl cellosolve (2-Butoxyethanol)	Tetrahydrofuran	Nitromethane
Butyl cellosolve acetate	Toluene	3-Pentanone (Diethyl ketone)
Carbon tetrachloride	1,1,1-Trichloroethane	Perchloroethylene
Cellosolve (2-Ethoxyethanol)	Propylene glycol butyl ether (PGBE)	a-Pinene
Cellosolve acetate	Ethyl lactate	n-Propanol
Chlorobenzene	Ethyl benzene	2-Propoxyethanol
Chloroform	1-Ethoxy-2-propyl acetate	n-Propyl acetate
o-Chlorotoluene	Heptane	Ethyl acrylate
Cumene	Hexane	1,1,2-Trichloroethane
Cyclohexanone	Isobutanol	Trichloroethylene
Diacetone alcohol	Isobutyl acetate	1,2,4-Trimethylbenzene
1,1-Dichloroethane	Isopropanol	1,3,5-Trimethylbenzene
1,2-Dichloroethane	Isopropyl acetate	Vinyl acetate
1,2-Dichloroethylene	d-Limonene	Xylene
Diisobutyl ketone	Methyl acetate	

PNAs (NIOSH 5506)

Acenaphthene	Benzo[b]fluoranthene	Chrysene	Indeno[1,2,3-cd]pyrene
Acenaphthylene	Benzo[k]fluoranthene	Dibenz[a,h]anthracene	Naphthalene
Anthracene	Benzo[ghi]perylene	Fluoranthene	Phenanthrene
Benz[a]anthracene	Benzo[a]pyrene	Fluorene	Pyrene

PNAs (OSHA 58)

Anthracene	Chrysene	Pyrene
Benzo[a]pyrene	Phenanthrene	

Sampling guide analyte descriptions and abbreviations

The following information is included in the Sample Guide's alphabetical listing of analytes:

Analyte Analytes are listed by their common name in alphabetical order in this sampling guide. The synonyms are listed within parenthesis. TVOC stands for Total Volatile Organic Compounds. VM&P Naphtha stands for Varnish Makers and Painters Naphtha.

CAS# Chemical Abstract Service number for the compound to be sampled.

Analytical method

This specifies the preferred analytical method used by Liberty Mutual for the analysis of the compound.

ASTM	American Society for Testing and Materials
Ciba-Geigy	Ciba-Geigy In-House Method for 1,3,5-triglycidyl isocyanurate
DuPont	DuPont In-House Method for Perfluorooctanoic Acid
HSE MDHS	UK Health & Safety Executive Methods for the Determination of Hazardous Substances
LM	Liberty Mutual Insurance "in-house" analytical methods
NIOSH	NIOSH Manual of Analytical Methods
OSHA	OSHA Manual of Analytical Methods

Analytical technique

AA	Atomic absorption spectrophotometry
AA-CV	Cold vapor atomic absorption spectrophotometry
EGA-TDA	Evolved gas analysis-thermo dilatometric analyzer
GC-ECD	Gas chromatography with electron capture detector
GC-FID	Gas chromatography with flame ionization detector
GC-MS	Gas chromatography-mass spectrometry
GC-NPD	Gas chromatography with nitrogen phosphorus detector
GC-TCD	Gas chromatography with thermal conductivity detector
GC-XSD	Gas chromatography with halogen specific detector
GFAA	Graphite furnace atomic absorption spectrophotometry
GRAV	Gravimetric analysis
HPLC	High performance liquid chromatography
IC	Ion chromatography
ICP	Inductively coupled plasma spectroscopy
ISE	Ion selective electrode
LC-MS	High performance liquid chromatography – mass spectrometry
PCM	Phase contrast microscopy
PLM	Polarized light microscopy
TEM	Transmission Electron Microscopy; Energy dispersive X-ray (EDX) analyzer
UV/VIS	UV/VIS spectrophotometry
XRD	X-ray diffractometry

Sampling media

The recommended sampling media for each of the methods in this sampling guide are:

AT571	Assay Technology Passive Monitor for Aldehydes
AT555	Assay Technology Passive Monitor for Ethylene Oxide
AT575	Assay Technology Passive Monitor for Nitrous Oxide
AT566	Assay Technology Passive Monitor for Organic Solvents
AgMF	25mm 0.45um Silver membrane filter
Anasorb708	SKC Anasorb 708 sorbent tube (SKC 226-30-08)
Anasorb747	SKC Anasorb 747 sorbent tube (SKC 226-81A)
Anasorb747/Anasorb747	SKC Anasorb 747 sorbent tube in series (SKC 226-82)
Bulk	Bulk sample
Cellulose Nitrate,Na₂CO₃	Cellulose nitrate filter treated with sodium carbonate (SKC 225-9031)
Carulite	Sorbent tube for mercury, replaces Hopcalite (SKC 226-17-1A/3A)
CS106	Chromosorb 106 sorbent tube (SKC 226-111A)
CT	Charcoal tube
CT, KOH	Potassium hydroxide treated Anasorb CSC coconut charcoal tube (SKC 226-67)
Cylinder	300cc Aluminum cylinder
di H₂O	Deionized water
GFF	Glass fiber filter
GFF, 1-2PP	1-(2-Pyridyl)piperazine treated glass fiber filter
GFF, Acid	Sulfuric acid treated glass fiber filter
GFF, HBr	Hydrogen bromide treated glass fiber filter
GFF, IOM	IOM sampler with glass fiber filter
GFF, NaNO₂	Sodium nitrite treated glass fiber filter
GFF,Vamine	Glass fiber filter coated with 10 mg of veratrylamine
GFF- PE,HgCl₂	Glass fiber filter -polyester filter treated with mercuric chloride in series (SKC 225-9018)
Ghost wipe	SKC wipe for surface lead or other metals
Impinger 4	Midget fritted glass bubbler containing 0.02% potassium iodide in sodium carbonate/sodium bicarbonate buffer
IOM	Personal inhalable sampler developed at the Institute of Occupational Medicine (IOM) in Scotland
MCE	Mixed cellulose ester filter membrane
MCE, 25mm 0.8um	Zefon PCM air sampling cassettes with conductive cowl (ZEFON Z008BA)
MCE,carbonate	Mixed cellulose ester filter treated with 20:1 sodium carbonate: glycerol solution
MCE-CS102	Mixed cellulose ester filter-Chromosorb 102 tube (SKC 226-104)

ORBO 34	Supelco specially treated charcoal tube for H ₂ S (20211)
ORBO 77	Supelco carbon beads treated with sulfuric acid for ammonia(20036)
ORBO 78	Supelco Carboxen-564 carbon molecular sieve treated with hydrogen bromide (20355)
ORBO 92	Supelco Carboxen-564 carbon molecular sieve (20362)
OVS-2/GFF	SKC XAD-2/glass fiber filter sorbent tube (226-30-16 OVS)
OVS-2/QF	SKC XAD-2/quartz filter sorbent tube (226-58 OVS)
OVS-7	SKC XAD-7/glass fiber filter sorbent tube (226-57 OVS)
OVS-Tenax	SKC tenax/glass fiber filter sorbent tube (226-56 OVS)
Oxidizer	SKC oxidizer with TEA-IMS for nitric oxide sampling
PS	SKC 520 inorganic mercury passive sampler (520-02A/03)
Paint chips	Paint chips sample in polyethylene bag
Porapak-P	SKC Porapak-P tube (226-114)
Porapak-Q	SKC Porapak-Q tube (226-115)
PTFE	25mm or 37mm 1.0um Polytetrafluoroethylene (PTFE) filter
PTFE-SGT	Teflon filter-silica gel tube in series
PTFE/XAD2	Teflon filter-XAD-2 tube in series
PTFE-AgMF	SKC polytetrafluoroethylene (Teflon) filter-silver membrane filter in series
PVC	37mm 5.0um pre-weighed polyvinyl chloride filter
PVC-KOH	Polyvinyl chloride filter- Midget fritted glass bubbler containing 0.1 N potassium hydroxide in series
QFF	Quartz fiber filter
QFF, titanium oxysulfate	Quartz fiber filter
QFF, Na₂CO₃	Quartz fiber filter treated with sodium carbonate
Sep-Pak	Waters dinitrophenylhydrazine (DNPH) treated cartridge for aldehydes (WAT047205)
SGT	Silica gel tube (SKC 226-10) with sorbent 75/150 mg
SGT**	SKC specially cleaned silica gel tube (226-10-03)
SGT/GFF	SKC silver nitrate coated silica gel tube with GFF coated with Na ₂ CO ₃ /glycerol) (SKC 226-177)
SGT,Acid	SKC sulfuric acid treated silica gel tube (226-10-06)
SGT, DNPH	SKC dinitrophenylhydrazine impregnated silica gel tube (226-119)
SGT***	2 specially washed and dried silica gel tubes in series (SKC 226-183)
SGT, HgCl₂	SKC mercuric chloride coated silica gel tube (226-10-02)
Soda Lime	SKC soda lime tube (226-210)
TEA-IMS	Triethanolamine-impregnated molecular sieve tube (SKC 226-40-02 / 226-40A with oxidizer)
Tenax	SKC tenax tube (226-35/226-35-01)
Whatman 42	Whatman 1442-070

XAD-2	SKC sorbent tube (226-30/226-30-04/226-30-06)
XAD-2,NITC	SKC sorbent tube coated with 10% 1-naphthylisothiocyanate (NITC) (226-30-18)
XAD-2,p-An	SKC sorbent tube (50/100) coated with 0.5mg of p-Anisidine (226-30-07)
XAD-7	SKC sorbent tube (226-95)
XAD-7, Acid	SKC XAD-7 sorbent (1,2-Dichloroethane) tube treated with 10% phosphoric acid (226-98)

Sampling rate (flow rate)

For passive monitors: cc/min. For all other compounds listed in the guide: liters per minute (lpm).

The sampling time in minutes for passive monitors and the sampling volume range indicates the minimum and maximum volume in liters. For bulk samples, the amount needed is listed in grams in this column.

Minimum volumes are typically calculated to allow quantification at 10% of the occupational exposure limits unless otherwise stated.

LOQ (Limit of Quantitation)

Limit of Quantitation is reported to two significant figures: %-Percent for bulk sample, fib/fld- Fibers per field, and µg- Micrograms/sample.

Compatibility code

This code indicates analytes that can be collected and analyzed simultaneously on a single sample. The absence of a code means that the analyte is not compatible with other analytes and would need to be collected on separate samples. The codes are as follows:

1%DMF/CS₂	Desorption in 1% dimethylformamide in carbon disulfide
1%IPA/CS₂	Desorption in 1% isopropanol in carbon disulfide
1%PRO/CS₂	Desorption in 1% n-propanol in carbon disulfide
5%IPA	Desorption in 5% isopropanol in deionized water
5% IPA/ CS₂	Desorption in 5% isopropanol in carbon disulfide
5%PRO/CS₂	Desorption in 5% n-propanol in carbon disulfide
95%EtOH	Desorption in 95% ethanol in deionized water
AC/CS₂	Desorption in 2% acetone in carbon disulfide
AC/MeOH	Desorption in 1% methanol in acetone
Acid1	Inorganic acids group 1 (see List of Scans for Inorganic Acids)
Acid2	Organic acids group 2 (formic acid, acetic acid, butyric acid, and propionic acid)
Acid3	Inorganic acids group on specialty media
Acid4	Inorganic acids group on specialty media
ACN	Desorption in acetonitrile
ACN/TOL	Desorption in 50% acetonitrile/toluene
ACN/DMSO	Desorption in 90% acetonitrile/dimethylsulfoxide
Acetone	Desorption in acetone
Aldehyde	Aldehyde group
Amine1	Aliphatic amine group by GC-FID
Amine2	Aliphatic amine group by HPLC

Amine3	Aromatic amine group
Benzene	Desorption in benzene
BUT/CS₂	Desorption in 1% 2-butanol in carbon disulfide
CCl₄	Desorption in carbon tetrachloride
Cl₂&Br₂	Chlorine and bromine
CS₂	Desorption in carbon disulfide
dil acid	Diluted sulfuric acid
DMF/CS₂	Desorption in 50% dimethylformamide in carbon disulfide
EA	Ethanolamine/diethanolamine/triethanolamine
Ethyl Acetate	Desorption in ethyl acetate
Ethyl Ether	Desorption in ethyl ether
FA	Desorption in formic acid
Isocyanate	Isocyanate and diisocyanate group
MC	Desorption in methylene chloride
MeOH	Desorption in methanol
MeOH/CS₂	Desorption in 1% methanol in carbon disulfide
%MeOH/MC	Desorption in mixture of methanol and methylene chloride
Metals	Metal group by either NIOSH 7301 or OSHA ID-125G
Metal1	Compatible metal group by modified NIOSH 7300 method
Metal2	Compatible metal group by modified OSHA ID-121 method
NO&NO₂	Nitric oxide and nitrogen dioxide
Phenol&cresol	Phenol and cresol group
PNAs	Polynuclear aromatic hydrocarbons
Silica	Silica by XRD
Tol	Desorption in toluene

Interferences

This column lists the possible interferences as stated in the analytical method.

Comments

This column contains special instructions for sample collection and handling.

Notes for gas analysis in comments

Note 1, Breathing air: Grade D and Grade E, Collected in Cylinders

Use when sampling a compressed gas line at 50 psig (i.e., before it goes into the regulator of the SCBA). If the pressure is lower than 50 psig, be sure to record the line pressure. This method reports CO, CO₂, %O₂, total hydrocarbons as CH₄ (TH), total halogenated hydrocarbons as CCl₄ (HH) and dew point (DP) in ppm and degrees F at 0 psig. See Procedure for Removing the Cylinder and LOQ at end of notes. Gravimetric analysis for condensed oil (total particulate) requires a 47mm Teflon filter.

Note 2, Medical gases – Contamination (NFPA 99)

Use when sampling a compressed medical gas line at 50 psig. When testing medical air, this method reports CO, CO₂, %O₂, total hydrocarbons as CH₄ (TH), total halogenated hydrocarbons as CCl₄ (HH) and dew point (DP) in ppm and degrees C at 50 psig. When testing nitrogen, this method reports CO, CO₂, total hydrocarbons as CH₄ (TH), total halogenated hydrocarbons as CCl₄ (HH) and dew point (DP) in ppm and degrees C at 50 psig. For oxygen, this method reports the same as medical air, except dew point and Total Hydrocarbons (TH). For nitrous oxide, this method reports CO, % Air and total halogenated hydrocarbons as CCl₄. See Procedure for Removing the Cylinder and LOQ at end of notes. Gravimetric analysis for condensed oil (total particulate) requires a 47-mm Teflon filter. The minimum sampling time for the gravimetric sample is 10 minutes in order to sample at least 1000 liters of gas, at 50 psig.

Procedure for removing the cylinder

A common mistake is removing the cylinder and the “A” fixture together. The correct steps are: Remove fixture “B.” The pop-off valve will activate. Remove the cylinder from fixture “A”, then remove fixture “A” from the gas line. Please refer to sampling instructions. Failure to follow the sampling instructions will result in the cylinder not being pressurized. ***Samples that are not pressurized cannot be analyzed.***

Limits of quantification

	Grav.	Air	CO	CO ₂	Dew Pt	O ₂	TH	HH
GAS	µg/s	%	ppm	ppm	see notes	%	ppm	ppm
Grade D&E	50	NA	0.50	25	200ppm, -33°F	0.30	0.50	0.50
Medical Air	50	NA	0.50	25	200ppm, -22°C	0.30	0.50	0.50
Nitrogen	50	NA	0.50	25	200ppm, -22°C	NA	0.50	0.50
Oxygen	50	NA	0.50	25	NA	0.30	N/A	0.50
Nitrous Oxide	50	0.15	0.50	NA	NA	NA	NA	0.50

Alphabetical list of analytes

18	Acenaphthene (see PNA scan)	23	Allyl Chloride
18	Acenaphthylene (see PNA scan)	23	Allyl Chloride
18	Acetaldehyde	23	Aluminum Metal and insoluble compounds
18	Acetaldehyde	23	Aluminum Oxide
19	Acetaldehyde	24	Ammonia
19	Acetic Acid	24	Ammonium Chloride Fume
19	Acetic Acid	24	Amyl Acetate
19	Acetone	24	Amyl Acetate
20	Acetone	25	Aniline
20	Acetonitrile	25	Anthracene (see PNA scan)
20	Acetonitrile	25	Anthracene (see PNA scan)
20	Acrolein	25	Antimony and compounds as Sb
21	Acrolein	26	Aromatic 100
21	Acrylamide	26	Arsenic and inorganic compounds, as As
21	Acrylic Acid	26	Arsenic and inorganic compounds, as As
21	Acrylonitrile (Vinyl Cyanide)	26	Arsine
22	Acrylonitrile (Vinyl Cyanide)	27	Asbestos, all forms
22	Aldehyde scan	27	Asbestos (bulk), all forms
22	Allyl Alcohol	27	Asbestos (Fiber Count)
22	Allyl Alcohol	27	Asphalt Fume
22	Allyl Bromide	28	Asphalt Fume as Benzene-Soluble Aerosol

28	Barium and soluble compounds as Ba	39	Butyl(n-) Alcohol
28	Benz[a]anthracene (see PNA scan)	40	Butyl(n-) Glycidyl Ether
28	Benzaldehyde	40	Butyl(n-) Glycidyl Ether
29	Benzaldehyde	40	Butyl(sec-) Acetate
29	Benzaldehyde	40	Butyl(sec-) Acetate
29	Benzene	40	Butyl(sec-) Alcohol
29	Benzene	41	Butyl(sec-) Alcohol
30	Benzo[a]pyrene (see PNA scan)	41	Butyl(tert-) Acetate
30	Benzo[a]pyrene (see PNA scan)	41	Butyl(tert-) Acetate
30	Benzo[b]fluoranthene (see PNA scan)	41	Butyl(tert-) Alcohol
30	Benzo[ghi]perylene (see PNA scan)	41	Butyl(tert-) Alcohol
31	Benzo[k]fluoranthene (see PNA scan)	42	Butyraldehyde(n-)
31	Benzyl Alcohol	42	Butyraldehyde(n-)
31	Benzyl Alcohol	42	Butyraldehyde(n-)
31	Benzyl Chloride	42	Butyric Acid
32	Benzyl Chloride	43	Cadmium and compounds as Cd
32	Beryllium and compounds as Be	43	Cadmium and compounds as Cd
32	Beryllium and Compounds as Be	43	Calcium
32	Biphenyl (Diphenyl)	44	Calcium Carbonate
33	Bismuth	44	Calcium Hydroxide
33	Bisphenol A	44	Calcium Oxide
33	Borate compounds, inorganic	45	Calcium Silicate Synthetic Nonfibrous
34	Breathing Air Grade D, Grade E	45	Calcium Sulfate (Gypsum)
34	Bromine	45	Camphor
34	Bromo(1-)-2-Chloroethane	45	Camphor
34	Bromoform	45	Caprolactam
35	Bromopropane(1-)	46	Carbaryl (SEVIN)
35	Bromopropane(1-)	46	Carbon Black
35	Butadiene(1,3-)	46	Carbon Black
35	Butadiene(1,3-)	46	Carbon Disulfide
36	Butanedione(2,3-); (Butadione(2,3-), Diacetyl, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)	47	Carbon Disulfide
36	Butanone(2-); (Methyl Ethyl Ketone)	47	Carbon Tetrachloride (tetrachloromethane)
36	Butanone(2-); (Methyl Ethyl Ketone)	47	Carbon Tetrachloride (tetrachloromethane)
36	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	48	Cellosolve (2-Ethoxyethanol)
37	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	48	Cellosolve (2-Ethoxyethanol)
37	Butoxyethoxy(2-(2-)) Ethanol	48	Ceramic Fibers (Fiber Count)
37	Butoxyethoxy(2-(2-)) Ethyl Acetate	48	Chlorine
37	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)	48	Chlorine Dioxide
37	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)	49	Chlorobenzene
38	Butyl Acrylate	49	Chlorobenzene
38	Butyl Acrylate	49	Chlorodiphenyl (Polychlorobiphenyl, 42% Chlorine)
38	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)	49	Chlorodiphenyl (Polychlorobiphenyl, 54% Chlorine)
39	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)	50	Chloroform (Trichloromethane)
38	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)	50	Chloroform (Trichloromethane)
38	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)	50	Chlorophenol(p-)
39	Butyl(n-) Acetate	50	Chloroprene(beta-); (2-Chloro- 1,3-butadiene)
39	Butyl(n-) Acetate	50	Chloroprene(beta-); (2-Chloro- 1,3-butadiene)
39	Butyl(n-) Alcohol	51	Chlorotoluene(o-)
		51	Chromium and Inorganic compounds as Cr

51	Chromium and Inorganic Compounds as Cr	62	Diesel Particulate Matter (Underground Mines)
52	Chromium, Hexavalent compounds as Cr	62	Diethanolamine
52	Chromium, Hexavalent Compounds as Cr	63	Diethylamine
52	Chrysene (see PNA scan)	63	Diethylenetriamine
53	Chrysene (see PNA scan)	62	Diethyl Ketone
53	Coal Dust – Anthracite	62	Diethyl Ketone (3- Pentanone)
53	Coal Dust – Bituminous	63	Diethyl Phthalate
53	Coal Tar Pitch Volatiles, as Benzene Soluble Aerosol	63	Diethyl Sulfate
54	Cobalt and Inorganic compounds as Co	63	Diglycidyl Ether of Bisphenol A
54	Cobalt and Inorganic compounds as Co	64	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)
54	Copper (Fume, Dusts and Mists) as Cu	65	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)
54	Copper (Fume, Dusts and Mists) as Cu	64	Dimethyl Acetamide
55	Cotton Dust	64	Dimethyl Acetamide
55	Cresol, all Isomers	64	Dimethyl Disulfide
55	Cumene	65	Dimethylformamide
55	Cumene	65	Dimethylformamide
56	Cyclohexane	64	Dimethyl Sulfide
56	Cyclohexane	65	Dioctyl Phthalate
56	Cyclohexanol	65	Dioxane(p-)
56	Cyclohexanol	66	Dioxane(p-)
56	Cyclohexanone	66	Diphenyl (Biphenyl)
57	Cyclohexanone	66	Dipropylene Glycol Methyl Ether Acetate (DPGMEA)
57	Cyclopentane	66	Dipropylene Glycol Methyl Ether (DPGME)
57	Desflurane (Suprane)	66	Dipropylene Glycol Methyl Ether (DPGME)
57	Desflurane (Suprane)	67	Divinyl Benzene
57	Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)	67	Enflurane (Ethrane)
58	Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)	67	Enflurane (Ethrane)
58	Diacetyl (Biacetyl, 2,3-Butadione, 2,3-Butanedione, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)	67	Epichlorohydrin (1-Chloro-2,3-epoxy propane)
58	Dibenzo[a,h]anthracene (see PNA scan)	67	Epichlorohydrin (1-Chloro-2,3-epoxy propane)
58	Dibutyl Ether	68	Ethanolamine (2-Aminoethanol)
59	Dibutyl Phthalate	68	Ethoxyethanol(2-) (Cellosolve)
59	Dichlorobenzene(o-)	68	Ethoxyethanol(2-) (Cellosolve)
59	Dichlorobenzene(o-)	68	Ethoxyethyl(2-) Acetate (Cellosolve acetate)
59	Dichlorobenzene(p-)	68	Ethyl 2-cyanoacrylate
59	Dichlorobenzene(p-)	69	Ethyl 3-ethoxypropionate
60	Dichloroethane(1,1)	69	Ethyl 3-ethoxypropionate
60	Dichloroethane(1,1)	69	Ethyl Acetate
60	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	69	Ethyl Acetate
60	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	69	Ethyl Acrylate
60	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)	70	Ethyl Acrylate
61	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)	70	Ethyl Alcohol (Ethanol)
61	Dichloromethane (Methylene chloride)	70	Ethyl Alcohol (Ethanol)
61	Dichloromethane (Methylene chloride)	71	Ethylamine
61	Diesel Particulate Matter (Above Ground Operations)	70	Ethyl Benzene
		70	Ethyl Benzene
		72	Ethylene Chlorohydrin (2-Chloroethanol)
		72	Ethylene Chlorohydrin (2-Chloroethanol)
		73	Ethylenediamine

72	Ethylene Dichloride (1,2-Dichloroethane)	84	Hydrogen Fluoride, as F
72	Ethylene Dichloride (1,2-Dichloroethane)	84	Hydrogen Peroxide
72	Ethylene Glycol	84	Hydrogen Sulfide
73	Ethylene Oxide	84	Hydroquinone (Dihydroxybenzene)
73	Ethylene Oxide	85	Hydroquinone (Dihydroxybenzene)
71	Ethyl Ether	85	Indium and Compounds as In
71	Ethyl Ether	85	Inorganic Acid Scan
71	Ethyl Lactate	85	Iodine and Iodides as I
71	Ethyl Methacrylate	86	Iron
73	Flour Dust	86	Iron Oxide
74	Fluoranthene (see PNA scan)	86	Isobutyl Acetate
74	Fluorene (see PNA scan)	86	Isobutyl Acetate
74	Fluorides, Particulate/Hydrogen Fluoride	87	Isobutyl Alcohol
74	Forane (Isoflurane)	87	Isobutyl Alcohol
75	Forane (Isoflurane)	87	Isocyanate Scan
75	Formaldehyde	87	Isoflurane (Forane)
75	Formaldehyde	88	Isoflurane (Forane)
75	Formaldehyde	88	Isooctane
76	Formamide	88	Isooctane
76	Formic Acid	88	Isophorone
76	Furfural	89	Isophorone
76	Furfuryl Alcohol	89	Isophorone Diisocyanate (IPDI)
77	Gasoline	89	Isophorone Diisocyanate (IPDI)
77	Gasoline	89	Isopropyl Acetate
77	Germanium	90	Isopropyl Acetate
77	Glutaraldehyde	90	Isopropyl Alcohol (Isopropanol)
78	Glutaraldehyde	90	Isopropyl Alcohol (Isopropanol)
78	Glutaraldehyde	90	Kaolin
78	Gold	91	Kerosene
78	Grain Dust	91	Kerosene
79	Graphite	91	Lactic Acid
79	Halothane (Fluothane)	91	Lanthanum
79	Halothane (Fluothane)	92	Lead and Inorganic Compounds as Pb
79	Heptane	92	Lead and Inorganic Compounds as Pb
80	Heptane	92	Lead and Inorganic Compounds as Pb
80	Heptanone(2-) (Methyl Amyl Ketone)	93	Lead Chromate as Cr(VI)
80	Heptanone(2-) (Methyl Amyl Ketone)	93	Limonene(d-)
81	Hexamethylene Diisocyanate (1,6-) (HDI)	93	Limonene(d-)
81	Hexamethylene Diisocyanate (1,6-) (HDI)	93	Lithium Salts
80	Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)	94	Magnesium
81	Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)	94	Magnesium Oxide
81	Hexane(n-)	94	Maleic Anhydride
82	Hexane(n-)	95	Manganese, elemental and Inorganic compounds as Mn
82	Hexyl Acrylate	95	Manganese, Elemental and Inorganic compounds as Mn
82	Hexylene Glycol (2-Methyl-2,4-pentanediol)	95	Medical Gases
82	Hydrazine	96	Mercury as Hg (Elemental and inorganic forms)
82	Hydrogen Bromide	96	Mercury as Hg (Elemental and inorganic forms)
83	Hydrogen Bromide	96	Mercury as Hg Particulate
83	Hydrogen Chloride	96	Metalworking Fluids
83	Hydrogen Chloride	97	
83	Hydrogen Cyanide		

97	Methanol (Methyl alcohol)	104	Methyl Methacrylate
97	Methanol (Methyl Alcohol)	104	Methyl Propyl Ketone (2-Pentanone)
97	Methanol (Methyl Alcohol)	104	Methyl Propyl Ketone (2-Pentanone)
98	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	104	Methyl Styrene(a-)
98	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	105	Methyl Styrene(a-)
98	Methoxyethanol(2-) (Methyl cellosolve, EGME)	105	Methyl Tert-butyl ether (MTBE)
98	Methoxyethanol(2-) (Methyl cellosolve, EGME)	105	Methyl Tert-butyl Ether (MTBE)
98	Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)	105	Methyl Vinyl Ketone
99	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	109	Mica
99	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	109	Mineral Oil (Oil mist)
99	Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)	109	Mineral Oil, used in metal working
106	Methyl(1-)-2-pyrrolidinone	109	Mineral Spirits (Stoddard Solvent)
106	Methyl(1-)-2-pyrrolidinone	110	Mineral Spirits (Stoddard Solvent)
99	Methyl Acetate	110	Molybdenum as Mo
99	Methyl Acetate	110	Molybdenum as Mo
100	Methyl Acrylate	110	Morpholine
100	Methyl Acrylate	111	Naphthalene
106	Methylacrylonitrile	111	Naphthalene
100	Methyl Alcohol (Methanol)	111	Naphthalene (see PNA scan)
100	Methyl Alcohol (Methanol)	111	Naproxen Sodium
101	Methyl Alcohol (Methanol)	112	Naproxen Sodium
101	Methyl Amyl Ketone (2-Heptanone)	112	Nickel and inorganic compounds as Ni
101	Methyl Amyl Ketone (2-Heptanone)	112	Nickel and inorganic compounds as Ni
101	Methyl Aniline	112	Nicotine
101	Methyl Chloroform (1,1,1-Trichloroethane)	113	Nitric acid
102	Methyl Chloroform (1,1,1-Trichloroethane)	113	Nitric Acid
106	Methylcyclohexane	113	Nitric Oxide and Nitrogen Dioxide
106	Methylcyclohexane	113	Nitroethane
102	Methyl Cyclopentane	114	Nitrogen Dioxide
108	Methylene(4,4')-bis(2 chloroaniline) (MOCA)	114	Nitromethane
108	Methylene(4,4') Dianiline (MDA)	114	Nitrous Oxide
107	Methylene Bis(4-cyclohexylisocyanate) (HMDI)	114	Oil Mist (Mineral Oil)
107	Methylene Bisphenyl Isocyanate (MDI)	115	Oil Mist (Mineral Oil) used in metal working
107	Methylene Bisphenyl Isocyanate (MDI)	115	Organic Solvent Scan
107	Methylene Bisphenyl Isocyanate, Polymeric (Polymeric MDI)	115	Oxalic Acid
108	Methylene Chloride (Dichloromethane)	115	Ozone
108	Methylene Chloride (Dichloromethane)	116	Palladium
102	Methyl Ethyl Ketone (2-Butanone, MEK)	116	Paraffin Wax Fume
102	Methyl Ethyl Ketone (2-Butanone, MEK)	116	Particles (insoluble or poorly soluble) Not otherwise specified; inhalable
103	Methyl Isoamyl Ketone	116	Particles (insoluble or poorly soluble) Not otherwise specified; respirable
103	Methyl Isoamyl Ketone	117	Particles (insoluble or poorly soluble) Not otherwise specified; total
103	Methyl Isobutyl Ketone (MIBK)	117	Pentanedione(2,3-)
103	Methyl Isobutyl Ketone (MIBK)	118	Pentanedione(2,4-)
103	Methyl Isopropyl Ketone	117	Pentane(n-)
		117	Pentane(n-)
		118	Pentanone(2-) (Methyl propyl ketone)
		118	Pentanone(2-) (Methyl propyl ketone)
		118	Peracetic Acid
		119	Perchloroethylene (Tetrachloroethylene)
		119	Perchloroethylene (Tetrachloroethylene)

119	Petroleum Ether	130	Pyridine
119	Petroleum Ether	130	Resin Acids
120	Phenanthrene (see PNA scan)	131	Resorcinol
120	Phenanthrene (see PNA scan)	131	Rhodium as Rh
120	Phenol	131	Scan for Aldehydes
121	Phenylcyclohexene (4-)	131	Scan for Inorganic Acids
120	Phenylcyclohexene(4-)	132	Scan for Isocyanates
121	Phenylene(1,3-) diamine	132	Scan for Organic Solvents
121	Phosphine	132	Scan for PNAs (NIOSH 5506)
121	Phosphoric Acid	132	Scan for PNAs (OSHA 58)
122	Phosphoric Acid	133	Selenium and Compounds as Se
122	Phosphorus (elements)	133	Sevoflurane (Sevofrane)
122	Phthalic Anhydride	133	Sevoflurane (Sevofrane)
122	Piperazine	133	Silica Cristobalite
123	Platinum Metal and Soluble Salts as Pt	134	Silica Quartz
123	PNA Scan (NIOSH 5506)	134	Silver Metal and Soluble Compounds as Ag
123	PNA Scan (OSHA 58)	134	Soapstone
124	Polychlorobiphenyl (Chlorodiphenyl, 42% Chlorine) (PCB)	134	Soapstone
123	Polychlorobiphenyl (Chlorodiphenyl, 54% Chlorine) (PCB)	135	Sodium
124	Polymeric (Methylene Bisphenyl Isocyanate) (Polymeric MDI)	135	Sodium Hydroxide
124	Polyvinyl Chloride (PVC)	135	Starch
124	Portland Cement	135	Stoddard Solvent
125	Potassium Hydroxide	136	Stoddard Solvent
125	Propanol(n-)	136	Strontium
125	Propanol(n-)	136	Styrene (Vinyl benzene)
125	Propionaldehyde	136	Styrene (Vinyl benzene)
126	Propionaldehyde	137	Sulfur Dioxide
126	Propionaldehyde	137	Sulfuric Acid
126	Propionaldehyde	137	Sulfuric Acid
126	Propionic Acid	137	Sulfuric Acid
126	Propoxyethanol(2-) (Ethylene glycol monopropyl ether)	138	Synthetic Vitreous Fibers (Fiber Count)
127	Propyl Bromide	138	Talc
127	Propyl Bromide	138	Tantalum and Tantalum Oxide Dust as Ta
129	Propylene Glycol (1,2-Propanediol)	138	Tellurium and Compounds as Te
128	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	139	Tetrachloroethylene (Perchloroethylene)
128	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	139	Tetrachloroethylene (Perchloroethylene)
128	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	139	Tetrahydrofuran
128	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	139	Tetrahydrofuran [THF]
129	Propylene Oxide (1,2-Epoxypropane)	140	Thallium and Compounds, as Tl
129	Propylene Oxide (1,2-Epoxypropane)	140	Thiram
127	Propyl(n-) Acetate	140	Tin and Compounds as Sn
127	Propyl(n-) Acetate	141	Titanium
127	Propyl(n-) Alcohol	141	Titanium Dioxide
128	Propyl(n-) Alcohol	141	Toluene
129	Pyrene (see PNA scan)	141	Toluene
130	Pyrene (see PNA scan)	142	Toluene-2,4-diisocyanate (2,4-TDI)
130	Pyrethrum	142	Toluene-2,4-diisocyanate (2,4-TDI)
		142	Toluene-2,6-diisocyanate (2,6-TDI)
		142	Toluene-2,6-diisocyanate (2,6-TDI)
		143	Toluidine(o-)
		143	Tributyl Phosphate
		143	Trichloro(1,1,2-)-1,2,2-trifluoroethane
		143	Trichloro(1,1,2-)-1,2,2-trifluoroethane

143 Trichlorobenzene(1,2,4-)
 144 Trichloroethane(1,1,1-) (Methyl Chloroform)
 144 Trichloroethane(1,1,1-) (Methyl Chloroform)
 144 Trichloroethane(1,1,2-)
 144 Trichloroethane(1,1,2-)
 144 Trichloroethylene
 145 Trichloroethylene
 145 Triethanolamine
 145 Triethylamine
 145 Triethylenetetramine
 146 Triglycidyl Isocyanurate(1,3,5)
 146 Trimellitic Anhydride
 146 Trimethylbenzene(1,2,4-)
 146 Trimethylbenzene(1,2,4-)
 146 Trimethylbenzene(1,3,5-)
 147 Trimethylbenzene(1,3,5-)
 147 Tungsten and Compounds as W
 (in the absence of Cobalt)
 147 Tungsten, as W Soluble Compounds
 147 TVOC as n-Hexane
 148 TVOC as n-Hexane
 148 Valeraldehyde
 148 Valeraldehyde
 148 Valeraldehyde
 149 Vanadium
 149 Vanadium Pentoxide as V
 149 Vegetable Oil Mist
 150 Vinyl(1-)-2-pyrrolidinone
 149 Vinyl Acetate
 150 Vinyl Acetate
 150 Vinyl Chloride (Chloroethylene)
 150 Vinyl Chloride (Chloroethylene)
 151 Vinylidene Chloride (1,1-Dichloroethylene)
 151 VM & P Naphtha
 151 VM & P Naphtha
 151 Welding Fume Scan
 152 Welding Fumes, Total
 152 Wood Dust
 152 Wood Dust
 152 Xylene (Dimethyl benzene)
 153 Xylene (Dimethyl benzene)
 153 Yttrium and compounds, as Y
 153 Zinc
 153 Zinc
 154 Zinc Chloride Fume
 154 Zinc Oxide

Acenaphthene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
83-32-9	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.33 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Acenaphthylene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
208-96-8	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.45 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Acetaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-07-0	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.7		15-480		0.029 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Acetaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-07-0	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		1-15		0.058 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Sample at 1 lpm for STEL.		

Acetaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-07-0	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		10-100		0.15 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm.		

Acetic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-19-7	NIOSH 2011	IC	PTFE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.5		5-100		2.3 µg
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Lab preferred method. Do not sample with inorganic acids.		

Acetic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-19-7	OSHA ID-PV2119	IC	Charcoal Tube (small)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.2		5-100		2.3 µg
Interferences		Comments		

Acetone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-64-1	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		0.5-3		0.91 µg
Interferences		Comments		
		Store and ship cold overnight.		

Acetone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-64-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
12.9		15-120		1.4 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Acetonitrile

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-05-8	NIOSH 1606	GC-FID	Charcoal Tube (large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		3-25		1.0 µg
Interferences		Compatibility Code		
		15%MeOH/MC		
Interferences		Comments		
Samples containing greater than 15% methanol or other alcohols		Large charcoal tubes are required for analyte collection since breakthrough volume is lower compared with smaller charcoal tubes. If also sampling for aldehydes, use AT571 monitors or Sep-Paks. Do not use DNPH tubes as they may off-gas acetonitrile.		

Acetonitrile

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-05-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
13.9		15-120		1.5 µg
Interferences		Compatibility Code		
		DMF/CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature. If also sampling for aldehydes, use AT571 monitors or Sep-Paks. Do not use DNPH tubes as they may off-gas acetonitrile.		

Acrolein

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-02-8	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.3		15-480		0.027 µg
Interferences		Compatibility Code		
		Aldehyde		
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved		Ship and store at controlled room temperature.		

Acrolein

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-02-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	10-100	0.14 µg	Aldehyde
Interferences		Comments	
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm.	

Acrylamide

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-06-1	OSHA 21	GC-FID	GFF/Silica Gel Tube (13 mm GFF/SKC 226-10)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	120	0.62 µg	MeOH
Interferences		Comments	
		Sample separately from CS ₂ compatible solvents.	

Acrylic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-10-7	OSHA PV2005	HPLC	Anasorb 708 (SKC 226-30-08) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	24	0.33 µg	
Interferences		Comments	
		Sample with 2 Anasorb 708 tubes in series. Separate and cap tubes after sampling.	

Acrylonitrile (Vinyl Cyanide)

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-13-1	OSHA 37	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	3.5-20	0.68 µg	CS ₂
Interferences		Comments	
		Preferred for STEL sampling. Sample at 0.2 lpm.	

Acrylonitrile (Vinyl Cyanide)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-13-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
12.9		15-480		1.0 µg
Interferences		Compatibility Code		
		DMF/CS ₂		
Interferences		Comments		

Aldehyde scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 2016	HPLC		
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
Interferences		Compatibility Code		
		See List of Scans for individual aldehydes.		

Allyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-18-6	NIOSH 1402	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		5-10		0.72 µg
Interferences		Compatibility Code		
		5% IPA/CS ₂		
Interferences		Comments		
		Sample separately from CS ₂ compatible solvents.		

Allyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-18-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
12.1		180-480		1.1 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		

Allyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-95-6	OSHA 1000	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-1.0		16-100		1.6 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Allyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-05-1	NIOSH 1000	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-1.0	16-100	0.96 µg	CS ₂
Interferences		Comments	
		Preferred for STEL sampling. Sample at a flow rate of 1 lpm.	

Allyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-05-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
10.5	180-480	1.4 µg	CS ₂
Interferences		Comments	

Aluminum Metal and insoluble compounds

CAS #	Analytical Method	Analytical Technique	Sampling Media
7429-90-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	200-1000	5.0 µg	Metals
Interferences		Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Aluminum Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
1344-28-1	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 LPM, for area sampling up to 15 LPM.	

Ammonia

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-41-7	OSHA ID-188	IC	Orbo 77 Tube
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-0.5	8-100	4.7 µg	
Interferences		Comments	
Particulate ammonium salts will interfere.		Sample at a flow rate of 0.5 lpm for STEL.	

Ammonium Chloride Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media
12125-02-9	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm; for area sampling up to 15 lpm.	

Amyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media
628-63-7	NIOSH 1450	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	2.3 µg	CS ₂
Interferences		Comments	

Amyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media
628-63-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.03	15-480	3.5 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Aniline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
62-53-3	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.02-0.5		5-30		1.0 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

Anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-12-7	OSHA 58	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		0.25 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

Anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-12-7	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.25 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Antimony and compounds as Sb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-36-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		50-1000		0.51 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Aromatic 100

CAS #	Analytical Method	Analytical Technique	Sampling Media
64742-95-6	NIOSH 1550	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	1.3-20	0.55 µg	CS ₂
Interferences	Comments		
Other aromatic compounds.	Please send bulk sample. Ship bulk samples separately from air samples.		

Arsenic and inorganic compounds, as As

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-38-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	560-1000	0.56 µg	Metals
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Arsenic and inorganic compounds, as As

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-38-2	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
NA	NA	6.5 µg	Metals2
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Arsine

CAS #	Analytical Method	Analytical Technique	Sampling Media
7784-42-1	NIOSH 6001	ICP	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	10	0.052 µg	
Interferences	Comments		
Other forms of Arsenic compounds (aerosol and gases) are quantified as Arsine.	High moisture may limit collection efficiency. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Asbestos (bulk), all forms

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-21-4	NIOSH 9002	PLM	Bulk	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
N.A.		1-10 grams		1.0%
Interferences		Comments		
Other fibers with optical properties similar to the asbestos minerals may give positive interferences.		Send bulk samples for asbestos analysis in double bagged ziplock bags with labels and chain of custody form outside the bag. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Asbestos (Fiber Count)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
Various	NIOSH 7400	PCM	25mm MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-16		50-720		0.050 asb/field
Interferences		Comments		
All other fibers will cause a positive interference. Chain-like particles may appear fibrous leading to high count and high levels of non-fibrous dust particles may obscure fibers leading to low count.		Sample open faced. Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. When shipping your samples, do not pack with untreated polystyrene as it can lead to fiber loss from electrostatic effect. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Asbestos, all forms

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-21-4	NIOSH 7402	TEM	25mm MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-16		50-720		0.050 asb/field
Interferences		Comments		
Other amphiboles similar to the asbestos minerals may give positive interferences. High dust background interferes with fiber identification.		Sample open faced. Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. This method is designed for use with Method 7400 (phase contrast microscopy). This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Asphalt Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-42-4	NIOSH 0500	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-15		40-7200		50 µg
Interferences		Comments		
All types of dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm; for area sampling up to 15 lpm.		

Asphalt Fume as Benzene-Soluble Aerosol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-42-4	NIOSH 5042	GRAV	Pre-weighed PTFE, 37mm	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		400-960		100 µg
Interferences		Comments		
All substances soluble in benzene will interfere. Changes in humidity pre- and post weighing can affect accuracy.		For inhalable sampling, please contact Lab for IOM sampler 1 week before intended use. The availability of IOM samplers is limited. Rental charges for use of IOM sampler apply.		

Barium and soluble compounds as Ba

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-39-3	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		0.50 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Benz[a]anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
56-55-3	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.082 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-52-7	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
6.97		15-480		0.047 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
100-52-7	NIOSH 2016	HPLC	Sep-Pak (WAT047205)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	10-100	0.24 µg	Aldehyde
Interferences		Comments	
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm for STEL	

Benzaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
100-52-7	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	1-15	0.094 µg	Aldehyde
Interferences		Comments	
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use flow rate of 1 lpm for STEL	

Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media
71-43-2	NIOSH 1501	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	5-30	0.40 µg	CS ₂
Interferences		Comments	
		Sample at flow rate of 0.2 lpm for STEL.	

Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media
71-43-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.69	15-480	0.60 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Benzo[a]pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-32-8	OSHA 58	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		0.20 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Benzo[a]pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-32-8	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.20 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Benzo[b]fluoranthene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
205-99-2	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.19 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Benzo[ghi]perylene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
191-24-2	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.27 µg
Interferences		Comments		
Asphalt fumes will interfere		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Benzo[k]fluoranthene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
207-08-9	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.25 µg
Interferences			Compatibility Code	
Asphalt fumes will interfere.			PNAs	
Interferences			Comments	
Asphalt fumes will interfere.			After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.	

Benzyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-51-6	NIOSH 1402	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-10		0.72 µg
Interferences			Compatibility Code	
			5%IPA/CS ₂	
Interferences			Comments	

Benzyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-51-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.50		15-480		1.1 µg
Interferences			Compatibility Code	
Interferences			Comments	
			Ship and store at controlled room temperature.	

Benzyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-44-7	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.73 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	

Benzyl Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media
100-44-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.14	15-480	1.1 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Beryllium and compounds as Be

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-41-7	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	400-1000	0.010 µg	Metals
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	Indicate in your sample submission form if beryllium oxide is suspected to be present. Beryllium oxide requires analysis by OSHA ID-125G. Minimum air volume required at 1/2 of TLV is 400 L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Beryllium and Compounds as Be

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-41-7	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
NA	NA	0.082 µg	Metals2
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	Indicate in your sample submission form if beryllium oxide is suspected to be present. Beryllium oxide requires analysis by OSHA ID-125G. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Biphenyl (Diphenyl)

CAS #	Analytical Method	Analytical Technique	Sampling Media
92-52-4	OSHA PV2022	GC-FID	XAD-7 (SKC-226-95)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.20	20	0.46 µg	CS ₂
Interferences	Comments		

Bismuth

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-69-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		0.50 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Bisphenol A

CAS #	Analytical Method	Analytical Technique	Sampling Media	
80-05-7	OSHA 1018	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.0		240		0.32 µg
Interferences		Comments		
		Samples must be stored and shipped cold.		

Borate compounds, inorganic

CAS #	Analytical Method	Analytical Technique	Sampling Media	
Varies	NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25		varies
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		LOQ for boron is 1.1 µg/sample. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. OEL is as inhalable. Order IOM samplers a week before survey date. Rental charges for IOM samplers apply.		

Breathing Air Grade D, Grade E

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0500 LMI-A5	GRAV GC-FID GC-TCD GC-XSD	47mm Teflon filter Cylinder	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		40-1000 (gravimetric) 300cc/25 psig (Cylinder)		Note 1
Interferences		Comments		
At high levels argon interferes with oxygen and nitrous oxide interferes with carbon dioxide.		A particulate blank is required.		

Bromine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7726-95-6	NIOSH 6011	IC	Silver Membrane Filter (with PTFE Pre-filter)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.3-1.0		70-360		4.5 µg
Interferences		Comments		
Hydrobromic acid will interfere.		Order media one week ahead, media is prepared when ordered.		

Bromo(1-)-2-Chloroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-04-0	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		1.0 µg
Interferences		Comments		

Bromoform

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-25-2	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-50		1.1 µg
Interferences		Comments		

Bromopropane(1-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
106-94-5	OSHA PV2061	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	12	1.3 µg	CS ₂
Interferences		Comments	

Bromopropane(1-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
106-94-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.26	15-480	2.0 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Butadiene(1,3-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
106-99-0	NIOSH 1024	GC-MS	CT-CT (SKC 226-01)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	10-25	1.1 µg	MC
Interferences		Comments	
High humidity (>80%) or other hydrocarbons present at permissible levels decrease sampler's capacity.		Use two large charcoal tubes in series. Separate and cap tubes after sampling. Ship cold overnight. If unable to ship overnight, store cold and then ship cold the following day.	

Butadiene(1,3-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
106-99-0	LM-GC-68	GC-MS	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
12.8	15-480	1.7 µg	MC
Interferences		Comments	
		Ship and store at controlled room temperature.	

Butanedione(2,3-); (Butadione(2,3-), Diacetyl, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)

CAS #	Analytical Method	Analytical Technique	Sampling Media
431-03-8	OSHA 1012	GC-ECD	Silica Gel Tube (SKC 226-183) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	9 (TWA) 3 (15-min short term)	0.040 µg	95% EtOH
Interferences	Comments		
	Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Separate and cap tubes after sampling. Order aluminum foil for wrapping the samples. For STEL sampling, a minimum of 3.0L is required. For TWA sampling, a minimum of 9.0L is required. Store and ship cold overnight. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Butanone(2-); (Methyl Ethyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
78-93-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
11.2	15-480	1.2 µg	CS ₂
Interferences	Comments		
Isopropyl acetate may co-elute with MEK.	Ship and store at controlled room temperature.		

Butanone(2-); (Methyl Ethyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
78-93-3	NIOSH 2500	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.25-12	0.79 µg	CS ₂
Interferences	Comments		
Isopropyl acetate may co-elute with MEK.	Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media
111-76-2	NIOSH 1403	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	2-10	1.2 µg	5%MeOH/MC
Interferences	Comments		

Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.43		15-480		1.8 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butoxyethoxy(2-(2-)) Ethanol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-34-5	OSHA PV2095	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		2-10		0.58 µg
Interferences		Compatibility Code		
		5%MeOH/MC		
Interferences		Comments		

Butoxyethoxy(2-(2-)) Ethyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
124-17-4	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.87 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.94 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.23		15-480		1.4 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-32-3	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.69 µg
Interferences		Comments		

Butyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-32-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.1		15-480		1.0 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.94 µg
Interferences		Comments		

Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-07-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.23		15-480		1.4 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl Cellosolve; (2-Butoxyethanol); (EGBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		2-10		1.2 µg
Interferences		Comments		
		5%MeOH/MC		

Butyl Cellosolve; (2-Butoxyethanol); (EGBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-76-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.43		15-480		1.8 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-86-4	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.84 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Butyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-86-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.94		15-480		1.3 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(n-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-36-3	NIOSH 1401	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		2-10		0.67 µg
Interferences		Compatibility Code		
		1%IPA/CS ₂		
Interferences		Comments		

Butyl(n-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-36-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.7		15-480		1.0 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(n-) Glycidyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2426-08-6	NIOSH 1616	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		15-30		1.1 µg
Interferences		Comments		

Butyl(n-) Glycidyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2426-08-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.7		15-480		1.7 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(sec-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
105-46-4	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.71 µg
Interferences		Comments		

Butyl(sec-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
105-46-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
		15-480		1.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(sec-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-92-2	NIOSH 1401	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		2-10		1.0 µg
Interferences		Comments		
		1%IPA/CS ₂		

Butyl(sec-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-92-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.7		15-480		1.5 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(tert-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-88-5	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.84 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Butyl(tert-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-88-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.50		15-480		1.3 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyl(tert-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-65-0	NIOSH 1400	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		2-10		0.55 µg
Interferences		Compatibility Code		
		BUT/CS ₂		
Interferences		Comments		

Butyl(tert-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-65-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.7		15-480		0.83 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Butyraldehyde(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.20		15-480	0.037 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Butyraldehyde(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5		10-100	0.19 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Sample at 1.0 lpm for STEL.		

Butyraldehyde(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-72-8	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5		1-15	0.073 µg	Aldehyde
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Sample at 1.0 lpm for STEL.		

Butyric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-92-6	NIOSH 2011	IC	PTFE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.5		15-100	2.3 µg	Acid2
Interferences		Comments		
		DO NOT sample with inorganic acids.		

Cadmium and compounds as Cd

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-43-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		500-1000		0.10 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Cadmium and compounds as Cd

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-43-9	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		0.50 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Calcium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-70-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		7.9 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Calcium Carbonate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1317-65-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		20 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.		

Calcium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1305-62-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		15 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.		

Calcium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1305-78-8	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		11 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of calcium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. If client wants different calcium salts speciated, please indicate in sample submission form or call/email Lab.		

Calcium Silicate Synthetic Nonfibrous

CAS #	Analytical Method	Analytical Technique	Sampling Media
1344-95-2	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences	Comments		
All other dusts will interfere.			

Calcium Sulfate (Gypsum)

CAS #	Analytical Method	Analytical Technique	Sampling Media
7778-18-9	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	20-7200	50 µg	
Interferences	Comments		
All other dusts will interfere.	For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

Camphor

CAS #	Analytical Method	Analytical Technique	Sampling Media
76-22-2	NIOSH 1301	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-25	0.53 µg	MeOH/CS ₂
Interferences	Comments		
	Preferred method. Sample at 0.2 lpm for STEL. Sample separately from CS ₂ compatible solvents.		

Camphor

CAS #	Analytical Method	Analytical Technique	Sampling Media
76-22-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
6.43	15-480	0.80 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Caprolactam

CAS #	Analytical Method	Analytical Technique	Sampling Media
105-60-2	OSHA PV2012	HPLC	OVS 7 (SKC 226-57)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	100	2.1 µg	
Interferences	Comments		

Carbaryl (SEVIN)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
63-25-2	OSHA 63	HPLC	OVS-2 (SKC 226-30-16)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		10-60		0.45 µg
Interferences			Comments	

Carbon Black

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1333-86-4	NIOSH 5000	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-2		30-570		50 µg
Interferences			Comments	
All other dusts will interfere.			Preferred method. For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.	

Carbon Black

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1333-86-4	OSHA ID-196	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		460-960		850 µg
Interferences			Comments	
Particulates that are insoluble in THF and that either vaporize or lose weight between 150°C and 600°C will interfere.			Please notify lab prior to sample collection.	

Carbon Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-15-0	NIOSH 1600	GC-MS	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		2-10		0.52 µg
Interferences			Comments	
			Sample separately from CS ₂ compatible solvents. Store and ship cold.	

Carbon Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-15-0	LM-GC-68	GC-MS	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.1		15-480		0.78 µg
Interferences		Comments		
		Sample separately from CS ₂ compatible solvents. Ship and store at controlled room temperature.		

Carbon Tetrachloride (tetrachloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
56-23-5	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.5-40		11 µg
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL. Preferred method for STEL sampling.		

Carbon Tetrachloride (tetrachloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
56-23-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.38		15-480		17 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Cellosolve (2-Ethoxyethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-6		0.54 µg
Interferences		Comments		
		5%MeOH/MC		

Cellosolve (2-Ethoxyethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-80-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.66		15-480		0.81 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ceramic Fibers (Fiber Count)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 7400	PCM	25mm MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-16		50-720		0.050 fiber/ field
Interferences		Compatibility Code		
Interferences		Comments		
All other fibers will cause a positive interference. Chain-like particles which may appear fibrous cause positive interference. High levels of non-fibrous dust particles may obscure fibers.		Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. Sample open faced. When shipping your samples, do not pack with untreated polystyrene as can lead to fiber loss from electrostatic effect. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Chlorine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7782-50-5	NIOSH 6011	IC	Silver Membrane Filter (with PTFE Pre-filter)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.3-1.0		80-360		2.3 µg
Interferences		Compatibility Code		
		Cl ₂ &Br ₂		
Interferences		Comments		
Hydrochloric acid will interfere.		Use a flow rate of 1.0 lpm for STEL. Order media one week ahead, media is prepared when ordered.		

Chlorine Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10049-04-4	OSHA ID-202	IC	Impinger 4	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5		35-120		0.75 µg
Interferences		Compatibility Code		
Interferences		Comments		
		Shelf-life of impinger solution is 3 months. Transfer sample solution to labeled vials before shipping. To avoid "freezing" of glass to glass, rinse impinger with distilled water before returning to the lab.		

Chlorobenzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-90-7	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.5-40		0.57 µg
Interferences			Comments	

Chlorobenzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-90-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.64		15-480		0.86 µg
Interferences			Comments	
			Ship and store at controlled room temperature.	

Chlorodiphenyl (Polychlorobiphenyl, 42% Chlorine)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
53469-21-9	NIOSH 5503	GC-MS	Florisil Tube & 13 mm GFF (SKC 226-39 with Swinnex Filter)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-50		0.97 µg
Interferences			Comments	
Other chlorinated pesticides may interfere in the quantification of PCB.			This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.	

Chlorodiphenyl (Polychlorobiphenyl, 54% Chlorine)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
11097-69-1	NIOSH 5503	GC-MS	Florisil Tube & 13 mm GFF (SKC 226-39 with Swinnex Filter)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-50		1.1 µg
Interferences			Comments	
Other chlorinated pesticides may interfere in the quantification of PCB.			This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.	

Chloroform (Trichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-66-3	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-50	5.9 µg	CS ₂
Interferences	Comments		

Chloroform (Trichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-66-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.39	15-480	8.9 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Chlorophenol(p-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
106-48-9	NIOSH 2014	HPLC	SGT (SKC 226-10)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	10-40	0.10 µg	
Interferences	Comments		

Chloroprene(beta-); (2-Chloro-1,3-butadiene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
126-99-8	NIOSH 1002	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.1	1.5-8	0.43 µg	CS ₂
Interferences	Comments		
	Store and ship cold.		

Chloroprene(beta-); (2-Chloro-1,3-butadiene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
126-99-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.34	15-480	0.65 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Chlorotoluene(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-49-8	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2		1.5-40	0.47 µg	CS ₂
Interferences		Comments		

Chromium and Inorganic Compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-47-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		670-1000	1.0 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of chromium are quantified. See chromium VI entry for special instructions on collection for this form of chromium. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Chromium and Inorganic compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-47-3	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
NA		NA	1.0 µg	Metals2
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of chromium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Chromium, Hexavalent compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-47-3	OSHA ID-215	IC	PVC or IOM-PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	500-960	0.050 µg	
Interferences		Comments	
Fe (II) appears to cause a negative interference during sampling and storage.		Please indicate in sample submission form the operation sampled (e.g., spray painting, chrome plating, welding, etc.) as the extraction method is different for paint samples. 500 liters is the minimum air volume at 50% TLV. 100 liters is the minimum air volume at 10% of the OSHA PEL. Ship all hexavalent chromium samples to the lab as soon as possible after sampling, via overnight shipping. OSHA method ID-215 states that samples from plating operations must be analyzed within 6 days of sampling and samples from welding operations must be analyzed within 8 days of sampling. If sampling CrVI as inhalable, order IOM samplers a week ahead of survey date.	

Chromium, Hexavalent Compounds as Cr

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-47-3	OSHA W4001	IC	QFF (Millipore AQFA03700) Wipes
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
NA	NA	0.050 µg	
Interferences		Comments	
Fe (II) appears to cause a negative interference during sampling and storage.		Please follow "Hexavalent Chromium Wipe Sampling Instructions." Do not use ghost wipes, Whatman, mixed cellulose ester (MCE) or glass fiber filters as they convert Cr(VI) to Cr(III). Please indicate in sample submission form the operation sampled as the extraction method is different for paint samples. Ship all hexavalent chromium samples to the lab as soon as possible after sampling, via overnight shipping. If sampling a chromium plating operation, request vials containing 10% Na ₂ CO ₃ with 2%NaHCO ₃ to stabilize the samples. Samples from welding operations must be analyzed within 8 days from date of sampling.	

Chrysene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media
218-01-9	OSHA 58	HPLC	Glass Fiber Filter
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	960	0.11 µg	PNAs
Interferences		Comments	
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.	

Chrysene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
218-01-9	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.11 µg
Interferences		Compatibility Code		
Asphalt fumes will interfere.		PNAs		
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Coal Dust – Anthracite

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.7		50-816		50 µg
Interferences		Compatibility Code		
All other respirable dusts will interfere.		PNAs		
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Coal Dust – Bituminous

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.7		50-816		50 µg
Interferences		Compatibility Code		
All other respirable dusts will interfere.		PNAs		
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Coal Tar Pitch Volatiles, as Benzene Soluble Aerosol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
65996-93-2	OSHA 58	GRAV	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		130 µg
Interferences		Compatibility Code		
The method is non-specific and measures all substances soluble in benzene.		PNAs		
Interferences		Comments		
The method is non-specific and measures all substances soluble in benzene.		This method is used to air monitor coke oven emissions, petroleum combustion products & asphalt fumes. If the BSF exceeds the PEL, then the sample is analyzed by HPLC to determine the presence of selected polynuclear aromatic hydrocarbons (PAHs). Wrap sample cassettes with aluminum foil after sampling. Ship and store cold.		

Cobalt and Inorganic compounds as Co

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-48-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		140-1000		0.10 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Cobalt and Inorganic compounds as Co

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-48-4	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		0.5 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Copper (Fume, Dusts and Mists) as Cu

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-50-8	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		70-1000		1.4 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Copper (Fume, Dusts and Mists) as Cu

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-50-8	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		5.0 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Cotton Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0500	GRAV	Pre-weighed PVC See Comment	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
See Comment		1000-3600		50 µg
Interferences		Comments		
All other dusts will interfere.		OSHA: Open-faced sampling with 3-piece pre-weighed PVC cassette on a vertical elutriator, cotton-dust sampler at 7.4 lpm. ACGIH: Use BGI GK2.69 cyclone at 1.6 lpm; sample 768 liters for 65% of new thoracic TLV.		

Cresol, all Isomers

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1319-77-3 95-48-7 108-39-4 106-44-5	OSHA 32	HPLC	XAD-7 (SKC 226-95)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		5-24		0.39 µg
Interferences		Comments		
		This method is applicable for all isomers of cresol (ortho-, meta-, and para-).		

Cumene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-82-8	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		2-30		0.55 µg
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.		

Cumene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-82-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.36		15-480		0.83 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Cyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-82-7	NIOSH 1500	GC-MS	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	2.5-5	0.17 µg	CS ₂
Interferences	Comments		

Cyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-82-7	LM-GC-68	GC-MS	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.82	15-360	0.25 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Cyclohexanol

CAS #	Analytical Method	Analytical Technique	Sampling Media
108-93-0	NIOSH 1402	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	1.1 µg	5%PRO/CS ₂
Interferences	Comments		

Cyclohexanol

CAS #	Analytical Method	Analytical Technique	Sampling Media
108-93-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.16	15-480	1.7 µg	MC
Interferences	Comments		
	Ship and store at controlled room temperature.		

Cyclohexanone

CAS #	Analytical Method	Analytical Technique	Sampling Media
108-94-1	NIOSH 1300	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	0.72 µg	CS ₂
Interferences	Comments		
	Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Cyclohexanone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-94-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.37		15-480		1.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
287-92-3	NIOSH 1500	GC-FID	CT (SKC 226-01)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		2.5-5.0		0.40 µg
Interferences		Comments		

Desflurane (Suprane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57041-67-5	OSHA 106	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05		3		6.1 µg
Interferences		Comments		
		Store and ship cold overnight.		

Desflurane (Suprane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57041-67-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.56		15-480		9.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-42-2	NIOSH 1402	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.98 µg
Interferences		Comments		
		5%PRO/CS ₂		

Diacetone Alcohol (4-Hydroxy-4-methyl-2-pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-42-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.5	15-480	1.5 µg	MC
Interferences	Comments		
	Ship and store at controlled room temperature.		

Diacetyl (Biacetyl, 2,3-Butadione, 2,3-Butanedione, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)

CAS #	Analytical Method	Analytical Technique	Sampling Media
431-03-8	OSHA 1012	GC-ECD	Silica Gel Tube (SKC 226-183) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	9 (TWA) 3 (15-min short term)	0.040 µg	95%EtOH
Interferences	Comments		
	Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Order aluminum foil for wrapping the samples during and after sampling. Separate and cap tubes after sampling. For STEL sampling, a minimum of 3.0L is required. For TWA sampling, a minimum of 9.0L is required. Store and ship cold overnight. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Dibenzo[a,h]anthracene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media
53-70-3	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	200-1000	0.93 µg	PNAs
Interferences	Comments		
Asphalt fumes will interfere.	After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Dibutyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media
142-96-1	NIOSH 1610	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	1-10	0.85 µg	CS ₂
Interferences	Comments		
	Store and ship cold.		

Dibutyl Phthalate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
84-74-2	NIOSH 5020	GC-FID	MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-3		10-200		0.46 µg
Interferences		Comments		
		Sample at 1.0 lpm for STEL.		

Dichlorobenzene(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-50-1	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.50 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Dichlorobenzene(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-50-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.55		15-480		0.75 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dichlorobenzene(p-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-46-7	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-8		0.45 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Dichlorobenzene(p-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-46-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.55		15-480		0.68 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dichloroethane(1,1)

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-34-3	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.5-15	1.7 µg	CS ₂
Interferences		Comments	
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.	

Dichloroethane(1,1)

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-34-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.22	15-480	2.5 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Dichloroethylene(1,2-)(trans); (Acetylene dichloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media
540-59-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.32	15-360	2.1 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Dichloroethylene(1,2-)(cis); (Acetylene dichloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media
540-59-0	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.2-5	1.6 µg	CS ₂
Interferences		Comments	

Dichloroethylene(1,2-)(cis); (Acetylene dichloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media
540-59-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.32	15-480	2.4 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Dichloroethylene(1,2-)(trans); (Acetylene dichloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-59-0	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		15-360		1.4 µg
Interferences		Comments		

Dichloromethane (Methylene chloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	NIOSH 1005	GC-FID	Charcoal tube (small) 2 Tubes in series	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		0.5-2.5		3.7 µg
Interferences		Comments		
		Separate tubes and cap immediately after sampling. Ship and store cold.		

Dichloromethane (Methylene chloride)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-09-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.54		15-240		5.6 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Diesel Particulate Matter (Above Ground Operations)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 5040	EGA-TDA	QFF	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2-4		142-19000		1.7 µg
Interferences		Comments		
		Order media at least one week ahead of time. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Diesel Particulate Matter (Underground Mines)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 5040	EGA-TDA	DPM cassettes (SKC 225-317)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7		142-19000	1.7 µg	
Interferences		Comments		
		Cyclone (SKC Part # 225-105) and cyclone holder (SKC Part # 225-1) should be used with these cassettes. Order media at least one week ahead of survey date. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Turnaround time is 10 business days.		

Diethanolamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-42-2	NIOSH 2007	IC	Silica Gel Tube (SKC 226-10-03)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.5		30	3.0 µg	EA
Interferences		Comments		
		Store in freezer after sampling. Ship cold.		

Diethyl Ketone (3- Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-22-0	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2		1-10	0.67 µg	CS ₂
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Diethyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-22-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.47		15-480	1.0 µg	CS ₂
Interferences		Comments		
		Ship and store at controlled room temperature.		

Diethyl Phthalate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
84-66-2	OSHA 104	GC-FID	OVS-Tenax (SKC 226-56)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		240		0.58 µg
Interferences		Comments		

Diethyl Sulfate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-67-5	NIOSH 2524	GC-FID	Porapak-P (SKC 226-114)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		0.25-12		7.2 µg
Interferences		Comments		

Diethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-89-7	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-1.0		5-30		0.31 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.		This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 15 business days.		

Diethylenetriamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-40-0	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		10		0.10 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

Diglycidyl Ether of Bisphenol A

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1675-54-3	OSHA 1018	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.5		270		0.50 µg
Interferences		Comments		

Dimethyl Acetamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-19-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.40		15-480		0.99 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dimethyl Acetamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-19-5	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-1		15-80		0.66 µg
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption.		

Dimethyl Disulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
624-92-0	LM-GC-59	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		5		1.6 µg
Interferences		Comments		
		Turnaround is 10 business days.		

Dimethyl Sulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-18-3	LM-GC-59	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		5		1.1 µg
Interferences		Comments		
		Turnaround is 10 business days.		

Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-83-8	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.59 µg
Interferences		Comments		

Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-83-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.69		15-480		0.89 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dimethylformamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
68-12-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.7		15-480		1.7 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dimethylformamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
68-12-2	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-1		15-80		1.5 µg
Interferences		Compatibility Code		
		MeOH		
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption.		

Diethyl Phthalate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
117-84-1	OSHA 104	GC-FID	OVS-Tenax (SKC 226-56)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.0		1-240		0.27 µg
Interferences		Compatibility Code		
		Tol		
Interferences		Comments		

Dioxane(p-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-91-1	NIOSH 1602	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-10		1.3 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Dioxane(p-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-91-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.77		15-480	2.0 µg	CS ₂
Interferences		Comments		
		Ship and store at controlled room temperature.		

Diphenyl (Biphenyl)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
92-52-4	OSHA PV2022	GC-FID	XAD-7 (SKC 226-95)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.20		20	0.46 µg	CS ₂
Interferences		Comments		

Dipropylene Glycol Methyl Ether (DPGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
34590-94-8	NIOSH 2554	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2		2-10	3.0 µg	MC
Interferences		Comments		

Dipropylene Glycol Methyl Ether (DPGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
34590-94-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
7.52		15-480	4.5 µg	MC
Interferences		Comments		
		Ship and store at controlled room temperature.		

Dipropylene Glycol Methyl Ether Acetate (DPGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
88917-22-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
6.37		15-480	1.5 µg	CS ₂
Interferences		Comments		
		Ship and store at controlled room temperature.		

Divinyl Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1321-74-0	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.57 µg
Interferences		Comments		

Enflurane (Ethrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13838-16-9	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-10		2.6 µg
Interferences		Comments		
		Store and ship cold overnight.		

Enflurane (Ethrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13838-16-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.15		15-480		4.0 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Epichlorohydrin (1-Chloro-2,3-epoxy propane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-89-8	NIOSH 1010	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		2-30		1.5 µg
Interferences		Comments		

Epichlorohydrin (1-Chloro-2,3-epoxy propane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-89-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.54		15-480		2.3 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethanolamine (2-Aminoethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
141-43-5	NIOSH 2007	IC	Silica Gel Tube (SKC 226-10-03)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.5	5	3.0 µg	EA
Interferences	Comments		
	Store in freezer after sampling. Ship cold.		

Ethoxyethanol(2-) (Cellosolve)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-80-5	NIOSH 1403	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	1-6	0.54 µg	5%MeOH/MC
Interferences	Comments		

Ethoxyethanol(2-) (Cellosolve)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-80-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.66	15-480	0.81 µg	MC
Interferences	Comments		
	Ship and store at controlled room temperature.		

Ethoxyethyl(2-) Acetate (Cellosolve acetate)

CAS #	Analytical Method	Analytical Technique	Sampling Media
111-15-9	NIOSH 1450	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	2.8 µg	CS ₂
Interferences	Comments		

Ethyl 2-cyanoacrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media
7085-85-0	OSHA 55	HPLC	XAD-7, Acid (SKC 226-98)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	12	0.70 µg	
Interferences	Comments		
	Ship and store cold.		

Ethyl 3-ethoxypropionate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
763-69-9	OSHA PV2025	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-10		1.1 µg
Interferences		Comments		

Ethyl 3-ethoxypropionate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
763-69-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.3		15-480		1.7 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-78-6	NIOSH 1457	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		1.2 µg
Interferences		Comments		

Ethyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
141-78-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.98		15-360		1.8 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
140-88-5	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.98 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. Store and ship cold overnight.		

Ethyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
140-88-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.16		15-480		1.5 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Alcohol (Ethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-17-5	NIOSH 1400	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		0.1-1		1.2 µg
Interferences		Comments		
		Store and ship cold overnight.		

Ethyl Alcohol (Ethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-17-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
13.5		15-120		1.8 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-41-4	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-30		0.5 µg
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.		

Ethyl Benzene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-41-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.91		15-480		0.75 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
60-29-7	NIOSH 1610	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.75 µg
Interferences		Comments		
Hexane may co-elute with the analyte of interest.		High humidity may greatly decrease breakthrough volume. Store and ship cold.		

Ethyl Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
60-29-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.3		15-240		1.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethyl Lactate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
687-47-8	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.85 µg
Interferences		Comments		

Ethyl Methacrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
97-63-2	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.65 µg
Interferences		Comments		
		Ship and store cold.		

Ethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-04-7	NIOSH 2010	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-1		5-30		NA
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.		This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 15 business days.		

Ethylene Chlorohydrin (2-Chloroethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-07-3	NIOSH 2513	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.51 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at flow rate of 0.2 lpm. High humidity may greatly decrease the breakthrough volume.		

Ethylene Chlorohydrin (2-Chloroethanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-07-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.2		15-480		0.77 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethylene Dichloride (1,2-Dichloroethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-06-2	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-50		1.5 µg
Interferences		Comments		

Ethylene Dichloride (1,2-Dichloroethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-06-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.22		15 -480		2.3 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Ethylene Glycol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-21-1	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-2		5-60		6.5 µg
Interferences		Comments		

Ethylene Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-2-8	ASTM D5578-04	GC-FID	AT555 Assay Technology Passive Monitor
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
12.4	15 – 480	0.69 µg	ACN/TOL
Interferences		Comments	
		Sample separately from CS ₂ compatible solvents. Ship and store at controlled room temperature.	

Ethylene Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-21-8	ASTM D5578-04	GC-FID	ORBO 78 (SUPELCO 20355)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.15	1-24	0.64 µg	ACN/TOL
Interferences		Comments	
		Store and ship cold. Sample separately from CS ₂ compatible solvents.	

Ethylenediamine

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-15-3	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	10	0.080 µg	Amine2
Interferences		Comments	
Nitrogen compounds that co-elute will interfere.			

Flour Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	HSE MDHS-14	GRAV	PVC, IOM
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	960	100 µg	
Interferences		Comments	
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab one week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.	

Fluoranthene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
206-44-0	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.30 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Fluorene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
86-73-7	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.32 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Fluorides, Particulate/Hydrogen Fluoride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 7906	IC	SKC 225-9031 (Cellulose Nitrate, Na ₂ CO ₃)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-2		70-1000		13 µg
Interferences		Comments		
Recovery of gaseous HF is reduced at high humidity.		Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.		

Forane (Isoflurane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
26675-46-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05		12		4.9 µg
Interferences		Comments		
		Store and ship cold overnight.		

Forane (Isoflurane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
26675-46-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.87		15-480		7.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-00-0	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
16.2		15-480		0.018 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-00-0	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		10-100		0.090 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use a flow rate of 1.0 lpm for STEL.		

Formaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-00-0	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		1-15		0.036 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Refrigerate media before and after sampling. Ship samples cold overnight. Preferred for STEL sampling. Use a flow rate of 1.0 lpm for STEL. For IAQ and LEED sampling where TVOC and formaldehyde are collected, do not sample with Charcoal Tubes and DNPH tubes in tandem. The sorbent in the DNPH tubes may off-gas acetonitrile which can cause a positive interference in the TVOC results.		

Formamide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-12-7	NIOSH 2004	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-1		15-80		3.1 µg
Interferences		Comments		
		Silica gel has a high affinity for water; high relative humidity may adversely affect the efficiency of analyte adsorption.		

Formic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
64-18-6	NIOSH 2011	IC	PTFE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.5		5-100		2.3 µg
Interferences		Comments		
		Do not sample with inorganic acids.		

Furfural

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-01-1	OSHA 72	GC-FID	SKC 226-81A	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1		180		0.96 µg
Interferences		Comments		
Furfuryl Alcohol is a sampling interference.				

Furfuryl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-00-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.87		15-480		0.96 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Gasoline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8006-61-9	NIOSH 1550	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.3-20		1.1 µg
Interferences		Comments		
		Please send bulk sample. Please ship bulk sample separately from the air samples. Preferred for STEL sampling. Use a flow rate of 0.2 lpm for STEL.		

Gasoline

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8006-61-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.19		15-480		1.7 µg
Interferences		Comments		
		Please send bulk sample. Please ship bulk sample separately from the air samples. Ship and store at controlled room temperature.		

Germanium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-56-4	NIOSH 7300	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		0.38 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-30-8	OSHA 64	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.24		15-480		0.011 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
111-30-8	NIOSH 2016	HPLC	Sep-Pak (WAT047205)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	10-100	0.055 µg	
Interferences	Comments		
	Refrigerate media before and after sampling. Ship cold overnight. Preferred for STEL sampling.		

Glutaraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
111-30-8	NIOSH 2532	HPLC	SGT, DNPH (SKC 226-119)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.5	1-30	0.022 µg	
Interferences	Comments		
	Refrigerate media before and after sampling. Ship cold overnight. Preferred for STEL sampling.		

Gold

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-57-5	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	50-1000	0.28 µg	Metals
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Grain Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	50-7200	50 µg	
Interferences	Comments		
All other dusts will interfere.	For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

Graphite

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7782-42-5	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.7		100-816		50 µg
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Halothane (Fluothane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
151-67-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05		12		3.4 µg
Interferences		Comments		
		Store and ship cold overnight.		

Halothane (Fluothane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
151-67-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.70		15-480		5.1 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Heptane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
142-82-5	NIOSH 1500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.45 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm.		

Heptane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
142-82-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.54		15-480		0.68 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Heptanone(2-) (Methyl Amyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	NIOSH 1301	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-25		0.64 µg
Interferences		Comments		

Heptanone(2-) (Methyl Amyl Ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-43-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.58		15-480		0.96 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28182-81-2	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		15		0.30 µg
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months from date prepared if kept cold. Ship cold overnight.		

Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28182-81-2	OSHA 5002	HPLC	Glass Fiber Filter Wipes	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
N/A		N/A	0.30 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.		

Hexamethylene Diisocyanate (1,6-) (HDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
822-06-0	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1		15	0.020 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months from date prepared if kept cold. Ship cold overnight.		

Hexamethylene Diisocyanate (1,6-) (HDI)

Hexamethylene Diisocyanate (HDI) (CAS# 822-06-0)				
CAS #	Analytical Method	Analytical Technique	Sampling Media	
822-06-0	OSHA 5002	HPLC	Glass Fiber Filter Wipes	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
N/A		N/A	0.020 µg	Isocyanate
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.		

Hexane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media		
110-54-3	NIOSH 1500	GC-FID	Charcoal Tube (small, large)		
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ	Compatibility Code
0.01-0.2		1-10		0.45 µg	CS ₂
Interferences			Comments		
			Preferred for STEL sampling. Sample at 0.2 lpm.		

Hexane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-54-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.88		15-480		0.72 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Hexyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2499-95-8	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.90 µg
Interferences		Comments		
		Store and ship cold overnight.		

Hexylene Glycol (2-Methyl-2,4-pentanediol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-41-5	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-2		5-60		2.5 µg
Interferences		Comments		

Hydrazine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
302-01-2	OSHA 108	HPLC	GFF, Acid	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		240		0.025 µg
Interferences		Comments		
		Media have short shelf-life. Media are prepared on request. Please contact the Lab 5 days before survey.		

Hydrogen Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10035-10-6	NIOSH 7903	IC	MCE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.2-0.5		10-100		4.6 µg
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Use a flow rate of 0.5 lpm for STEL sampling.		

Hydrogen Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media
10035-10-6	NIOSH 7907	IC	SKC 225-9032
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	35-600	3.8 µg	Acid3
Interferences		Comments	
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results.		Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.	

Hydrogen Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media
7647-01-0	NIOSH 7903	IC	MCE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-0.5	10-100	2.3 µg	Acid1
Interferences		Comments	
Particulate salts of the acid will give a positive interference.		Use a flow rate of 0.5 lpm for STEL sampling.	

Hydrogen Chloride

CAS #	Analytical Method	Analytical Technique	Sampling Media
7647-01-0	NIOSH 7907	IC	SKC 225-9032
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	40-600	13 µg	Acid3
Interferences		Comments	
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results. Potentially interfering particulate chlorides and nitrates removed by the pre-filter can react with the sampled acids and liberate HCl and HNO ₃ , which gets collected on the sampling filter, leading to high results.		Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.	

Hydrogen Cyanide

CAS #	Analytical Method	Analytical Technique	Sampling Media
74-90-8	NIOSH 7904	ISE	Soda Lime Tube (SKC 226-210)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2	10-90	5.2 µg	
Interferences		Comments	
Sulfide, chloride, iodide, bromide, cadmium, zinc, silver, nickel, cuprous iron and mercury interfere.		Use a flow rate of 0.2 lpm for STEL sampling. Method is not covered under our AIHA LAP, LLC scope of accreditation. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 15 business days.	

Hydrogen Fluoride, as F

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-39-3	NIOSH 7903	IC	MCE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-0.5	30-100	3.0 µg	Acid1
Interferences		Comments	
Particulate salts of the acid will give a positive interference.		Use the maximum flow rate for STEL sampling.	

Hydrogen Peroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
7722-84-1	OSHA ID-1019	UV/VIS	SKC 225-9030 (QFF, titanium oxysulfate)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.0	240	10.0 µg	
Interferences		Comments	
Any compound with a response, or that reacts with the titanium oxysulfate to produce a response, at 410nm is a potential interferent.		Refrigerate media before and after sampling. Ship samples cold overnight. After sampling seal the cassette with the end plugs and wrap each cassette with foil. Order foil for wrapping the samples. Filters have limited shelf-life. This method is not covered under our AIHA LAP, LLC scope of accreditation. Specialty Filter. Media charge applies.	

Hydrogen Sulfide

CAS #	Analytical Method	Analytical Technique	Sampling Media
7783-06-4	NIOSH 6013	IC	ORBO 34 (SUPELCO 20211)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	30-70	8.9 µg	
Interferences		Comments	
Sulfur dioxide gas may give a positive interference for hydrogen sulfide.		Use a flow rate of 0.5 lpm for STEL.	

Hydroquinone (Dihydroxybenzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-31-9	NIOSH 5004	HPLC	MCE Filter
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	30-180	10 µg	
Interferences		Comments	
		Hydroquinone is unstable on the collecting media. Stabilize immediately after collecting by transferring filter into a vial containing 1% acetic acid.	

Hydroquinone (Dihydroxybenzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-31-9	OSHA PV2094	HPLC	XAD-7, Acid (SKC 226-98)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2	20	0.30 µg	
Interferences	Comments		
	Preferred method.		

Indium and Compounds as In

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-74-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	150-1000	1.5 µg	Metals
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Inorganic Acid Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 7903	IC	MCE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-0.5	50-100	N/A	
Interferences	Comments		
	See List of Scans for list of individual inorganic acid/anion.		

Iodine and Iodides as I

CAS #	Analytical Method	Analytical Technique	Sampling Media
7553-56-2	NIOSH 6005	IC	CT, KOH (SKC 226-67)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.5-1.0	50-200	5.4 µg	
Interferences	Comments		
Particulate iodide salts, hydrogen iodide or organic iodides may give a positive interference.	Use a flow rate of 1.0 lpm for STEL.		

Iron Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-37-1	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PCV	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		4.3 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of iron are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Iron

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-37-1	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		30 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Isobutyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-19-0	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.81 µg
Interferences		Comments		

Isobutyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-19-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.5		15-480		1.2 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Isobutyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media
78-83-1	NIOSH 1401	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	0.59 µg	1%PRO/CS ₂
Interferences	Comments		

Isobutyl Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media
78-83-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
10.7	15-480	0.89 µg	MC
Interferences	Comments		
	Ship and store at controlled room temperature.		

Isocyanate Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	15		
Interferences	Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.	See List of Scans for individual isocyanates. Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months from date prepared if kept cold. Ship cold overnight.		

Isoflurane (Forane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
26675-46-7	OSHA 103	GC-FID	Anasorb 747 (SKC 226-81A)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05	12	4.9 µg	CS ₂
Interferences	Comments		
	Store and ship cold overnight.		

Isoflurane (Forane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
26675-46-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.87		15-480		7.4 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Isooctane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-84-1	NIOSH 1500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.52 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Isooctane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
540-84-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.60		15-480		0.78 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Isophorone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-59-1	NIOSH 2508	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.90 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL. High humidity may greatly decrease the breakthrough volume.		

Isophorone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-59-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.8		15-480		1.4 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Isophorone Diisocyanate (IPDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4098-71-9	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		15		0.034 µg
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Sample open-faced. Media must be stored cold before and after sampling. Filter is stable for 6 months from date prepared if kept cold. Ship cold overnight.		

Isophorone Diisocyanate (IPDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4098-71-9	OSHA 5002	HPLC	Glass Fiber Filter Wipes	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
N/A		N/A		0.034 µg
Interferences		Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.		

Isopropyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-21-4	NIOSH 1454	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.02-0.2		0.1-9		1.0 µg
Interferences		Comments		
		Use a flow rate of 0.2 lpm for STEL.		

Isopropyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-21-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.08		15-480		1.5 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Isopropyl Alcohol (Isopropanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-63-0	NIOSH 1400	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		0.3-3		1.0 µg
Interferences		Compatibility Code		
		BUT/CS ₂		
Interferences		Comments		

Isopropyl Alcohol (Isopropanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
67-63-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.8		15-240		1.5 µg
Interferences		Compatibility Code		
		ACN		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Kaolin

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1332-58-7	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.7		100-816		50 µg
Interferences		Compatibility Code		
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Kerosene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8008-20-6	NIOSH 1550	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.3-20		2.9 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

Kerosene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8008-20-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
6.42		15-480		4.4 µg
Interferences		Comments		
		Ship and store at controlled room temperature. Please send bulk sample. Ship bulk sample separately from air samples.		

Lactic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
50-21-5	NIOSH 2011	IC	PTFE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.2-0.5		5-100		3.0 µg
Interferences		Comments		
		Do not sample for inorganic acids using the same tube.		

Lanthanum

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-91-0	NIOSH 7301	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		120-1000		0.10 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-92-1	OSHA ID-121	ICP	ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		1.0 µg
Interferences			Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-92-1	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		120-1000		0.25 µg
Interferences			Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Lead and Inorganic Compounds as Pb

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-92-1	OSHA ID-121 NIOSH 7301	ICP	Paint chips	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		0.25 µg
Interferences			Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			Needs at least 1 gram of bulk sample. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Lead Chromate as Cr(VI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7758-97-6	OSHA ID-215	IC	PVC or IOM-PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		500-960		0.05 µg
Interferences		Comments		
		Please indicate in sample submission form the operation sampled (e.g., spray painting, chrome plating, welding, etc.) as the extraction method is different for paint samples. 500 liters is the minimum air volume at 50% TLV. 100 liters is the minimum air volume at 10% of the OSHA PEL. Ship all hexavalent chromium samples to the lab as soon as possible after sampling, via overnight shipping. OSHA method ID-215 states that samples from plating operations must be analyzed within 6 days of sampling and samples from welding operations must be analyzed within 8 days of sampling. If sampling CrVI as inhalable, order IOM samplers a week ahead of survey date.		

Limonene(d-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
138-86-3	NIOSH 1552	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.52 µg
Interferences		Comments		
		CS ₂		

Limonene(d-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
138-86-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.86		15-480		0.78 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Lithium Salts

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-93-1	NIOSH 7301	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		250-1000		0.10 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of lithium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Magnesium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-95-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		25-1000	1.0 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of magnesium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Magnesium Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1309-48-4	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		25-1000	1.7 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of magnesium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Maleic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-31-6	OSHA 86	HPLC	GFF, Vamine	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.5		60	0.74 µg	ACN/DMSO
Interferences		Comments		
Both phthalic amd trimellitic anhydride should be considered as potential sampling interferences.		Sampling media has short shelf-life, so media is prepared when ordered. Please order filters at least 5 days prior to survey date.		

Manganese, Elemental and Inorganic compounds as Mn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-96-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		125-1000		0.15 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. Manganese may be sampled using respirable PPI. Use a flow rate of 2 lpm. Please indicate in your request if you want the filter inside the PPI pre-weighed. Please contact lab a week before your survey if using respirable PPI. Media charge applies for PPIs.		

Manganese, elemental and Inorganic compounds as Mn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-96-5	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		0.60 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Medical Gases

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0500 LMI-A5	GRAV, GC-FID, GC-TCD GC-XSD	47mm Teflon filter Cylinder	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
Note 2* see Comments		1000 (gravimetric) 300cc/25 psig (Cylinder)		Note 2* see Comments
Interferences		Comments		
At high levels argon interferes with oxygen and carbon dioxide interferes with nitrous oxide.		*See Note 2 in the Comments section of "IH Lab Sampling Guide Analyte Descriptions and Abbreviations"		

Mercury as Hg (Elemental and inorganic forms)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	NIOSH 6009	AA-CV	Carulite Tube (SKC 226-17-1A/3A)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.15-0.25		50-120		0.12 µg
Interferences		Comments		
Inorganic and organic mercury compounds may cause a positive interference. Oxidizing gases, including chlorine, do not interfere.		A 37-mm, cellulose ester membrane filter in a cassette preceding the sorbent may be used if particulate mercury is to be determined separately. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Mercury as Hg (Elemental and inorganic forms)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	OSHA ID-140	AA-CV	PS (SKC 520-02A/03)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.020		9.6		0.12 µg
Interferences		Comments		
Particulate mercury compounds are a positive interference.		Call Lab one week before sampling. Lab has limited number of passive mercury monitor holders. Refer to OSHA ID-145 if sampling in workplaces containing both mercury vapor and particulate. Specialty media. Media charge applies. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Mercury as Hg Particulate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-97-6	OSHA ID-145	AA-CV	Air MCE (SKC225-5) Wipe (Whatman#42) Bulk	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		50		0.12 µg
Interferences		Comments		
		For wipe sampling, moisten the wipe filter with deionized water prior to use. If mercury vapor is suspected to be present, please refer to OSHA ID-140 or NIOSH 6009 for additional sampling information. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Metalworking Fluids

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 5524	Extraction/GRAV	Pre-weighed PTFE, 37mm
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	768-960	75 µg	
Interferences		Comments	
The method is non-specific and measures all substances extractable by organic solvents.		NIOSH 5524 recommends submitting one bulk sample of each type of fluid for solubility testing. After sampling, ship and store cold. Please note that samples must be analyzed within 2 weeks after collection.	

Methanol (Methyl alcohol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	Anasorb 747 (SKC 226-82) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05	1-5	2.3 µg	DMF/CS ₂
Interferences		Comments	
		Sample with 2 Anasorb 747 tubes in series ("Part A" as the front section and "Part B" as the back section. Please order as a set.) Recommended air volume is 5 liters when relative humidity is >50% and 3 liters when relative humidity is <50% at 25°C. Separate and cap tubes after sampling. Store and ship cold overnight.	

Methanol (Methyl Alcohol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	SKC 575-002 Passive Sampler
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
24	15-60	2.3 µg	DMF/CS ₂
Interferences		Comments	
		Use for STEL sampling. Store and ship cold overnight.	

Methanol (Methyl Alcohol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	SKC 575-007 Passive Sampler
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.2	60-480	2.3 µg	DMF/CS ₂
Interferences		Comments	
		Green Diffusive Barrier must be used when sampling. Store and ship cold overnight.	

Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-10		1.3 µg
Interferences		Comments		

Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.66		15-480		2.0 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methoxyethanol(2-) (Methyl cellosolve, EGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-86-4	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-50		1.7 µg
Interferences		Comments		
		Preferred method for STEL sampling. Sample at 0.05 lpm.		

Methoxyethanol(2-) (Methyl cellosolve, EGME)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-86-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
10.5		120-480		2.6 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
111-77-3	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-10		1.3 µg
Interferences		Comments		
		5%MeOH/MC		

Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-49-6	NIOSH 1451	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.2-20	0.46 µg	CS ₂
Interferences		Comments	
		Preferred method for STEL sampling.	

Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-49-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.00	120-480	0.67 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)

CAS #	Analytical Method	Analytical Technique	Sampling Media
111-96-6	NIOSH 1403	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	1-10	0.58 µg	5%MeOH/MC
Interferences		Comments	

Methyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-20-9	NIOSH 1458	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.2-10	1.5 µg	CS ₂
Interferences		Comments	
		Preferred method for STEL sampling. Sample at 0.2 lpm.	

Methyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-20-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
10.3	15-480	2.3 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Methyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media
96-33-3	NIOSH 1459	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-5	1.1 µg	CS ₂
Interferences		Comments	

Methyl Acrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media
96-33-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.89	15-480	1.7 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Methyl Alcohol (Methanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	Anasorb 747 (SKC 226-82) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05	1-5	2.3 µg	DMF/CS ₂
Interferences		Comments	
		Sample Anasorb 747 tubes in series (with "Part A" in front and "PART B" in back section. Recommended air volume is 5 liters when relative humidity is >50% and 3 liters when relative humidity is <50% at 25° C. Separate and cap tubes after sampling. Store and ship cold overnight.	

Methyl Alcohol (Methanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	SKC 575-002 Passive Sampler
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
24	15-60	2.3 µg	DMF/CS ₂
Interferences		Comments	
		Use for STEL sampling. Store and ship cold overnight.	

Methyl Alcohol (Methanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media
67-56-1	OSHA 5001	GC-FID	SKC 575-007 Passive Sampler
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.2	60-480	2.3 µg	DMF/CS ₂
Interferences	Comments		
	Green Diffusive Barrier must be used when sampling. Store and ship cold overnight.		

Methyl Amyl Ketone (2-Heptanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-43-0	NIOSH 1301	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-25	0.61 µg	MeOH/CS ₂
Interferences	Comments		
	Use a flow rate of 0.2 lpm for STEL.		

Methyl Amyl Ketone (2-Heptanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-43-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.58	15-480	0.92 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Methyl Aniline

CAS #	Analytical Method	Analytical Technique	Sampling Media
100-61-8	NIOSH 2002	GC-FID	SGT (SKC 226-10)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.02-0.5	5-30	NA	Amine3
Interferences	Comments		
Nitrogen compounds that co-elute will interfere.			

Methyl Chloroform (1,1,1-Trichloroethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
71-55-6	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	2.2 µg	CS ₂
Interferences	Comments		

Methyl Chloroform (1,1,1-Trichloroethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-55-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.93		15-480		3.3 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methyl Cyclopentane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
96-37-7	NIOSH 1500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.38 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Methyl Ethyl Ketone (2-Butanone, MEK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	NIOSH 2500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		0.25-12		0.86 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. Ship and store cold.		

Methyl Ethyl Ketone (2-Butanone, MEK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-93-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.2		15-480		1.3 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
Isopropyl acetate may co-elute with MEK.		Ship and store at controlled room temperature.		

Methyl Isoamyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-12-3	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.66 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	

Methyl Isoamyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-12-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.58		15-480		0.99 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	
			Ship and store at controlled room temperature.	

Methyl Isobutyl Ketone (MIBK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-10-1	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.54 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	
			Preferred for STEL sampling. Sample at 0.2 lpm.	

Methyl Isobutyl Ketone (MIBK)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-10-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.16		15-480		0.80 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	
			Ship and store at controlled room temperature.	

Methyl Isopropyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
563-80-4	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		1.0 µg
Interferences			Compatibility Code	
			CS ₂	
Interferences			Comments	
			Ship and store cold.	

Methyl Methacrylate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
80-62-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.16		15-480		1.5 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methyl Propyl Ketone (2-Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.73 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Methyl Propyl Ketone (2-Pentanone)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-87-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.89		15-480		1.1 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methyl Styrene(a-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-83-9	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-30		0.41 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. Store and ship cold overnight.		

Methyl Styrene(a-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
98-83-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.43		15-480		0.62 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methyl Tert-butyl Ether (MTBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1634-04-4	NIOSH 1615	GC-FID	Charcoal Tube, small (2 Tubes in series)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-0.2		2-96		0.67 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Preferred for STEL sampling. Sample at 0.2 lpm. Store and ship cold immediately.		

Methyl Tert-butyl ether (MTBE)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1634-04-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.77		15-480		1.0 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methyl Vinyl Ketone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
78-94-4	NIOSH 1300	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.27 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Store and ship cold overnight.		

Methyl(1-)-2-pyrrolidinone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
872-50-4	NIOSH 1302	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		0.5-125		1.1 µg
Interferences		Compatibility Code		
		5%MeOH/MC		
Interferences		Comments		

Methyl(1-)-2-pyrrolidinone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
872-50-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.31		15-480		1.7 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methylacrylonitrile

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-98-7	NIOSH 1604	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		3.5-20		1.0 µg
Interferences		Compatibility Code		
		AC/CS ₂		
Interferences		Comments		

Methylcyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-87-2	NIOSH 1500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		4		0.45 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Methylcyclohexane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-87-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.32		15-420		0.68 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Methylene Bis(4-cyclohexylisocyanate) (HMDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
5124-30-1	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		15		0.070 µg
Interferences			Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.			Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.	

Methylene Bisphenyl Isocyanate (MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-68-8	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		15		0.014 µg
Interferences			Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.			Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.	

Methylene Bisphenyl Isocyanate (MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
101-68-8	OSHA 5002	HPLC	Glass Fiber Filter Wipes	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
N/A		N/A		0.014 µg
Interferences			Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.			Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.	

Methylene Bisphenyl Isocyanate, Polymeric (Polymeric MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
9016-87-9	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		15		0.10 µg
Interferences			Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.			Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.	

Methylene Chloride (Dichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-09-2	NIOSH 1005	GC-FID	CT-CT (SKC 226-01)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	0.5-2.5	2.8 µg	CS ₂
Interferences		Comments	
		Sample using two CT in series. Separate and cap tubes after sampling. Ship and store cold immediately.	

Methylene Chloride (Dichloromethane)

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-09-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.54	15-240	4.2 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Methylene(4,4')-Dianiline (MDA)

CAS #	Analytical Method	Analytical Technique	Sampling Media
101-77-9	NIOSH 5029	HPLC	GFF, Acid
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-2	10-1000	0.68 µg	
Interferences		Comments	
4,4'-Diphenyl methane diisocyanate (MDI) will interfere.		Media has short shelf life so it is not kept in stock. Please contact the lab one week before sampling to order media. Within 4 hours of sampling, transfer the filter to a glass vial containing 4 ml 0.1 N methanolic KOH. This method can also be used for wipe sampling.	

Methylene(4,4')-bis(2 chloroaniline) (MOCA)

CAS #	Analytical Method	Analytical Technique	Sampling Media
101-14-4	NIOSH P&CAM 236	HPLC	GFF-SGT (SKC 225-16, SKC226-10)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-1	50	0.50 µg	
Interferences		Comments	

Mica

CAS #	Analytical Method	Analytical Technique	Sampling Media
12001-26-2	NIOSH 0600	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7	500-1200	50 µg	
Interferences	Comments		
All other respirable dusts will interfere.	Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and BGI-4L at 2.2 lpm. 2020 NIC of 0.1mg/m ³ (R) was adopted in 2021.		

Mineral Oil (Oil mist)

CAS #	Analytical Method	Analytical Technique	Sampling Media
8012-95-1	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences	Comments		
All other dusts will interfere.	This method is not for oil mist containing PNAs; it does not collect vapor.		

Mineral Oil, used in metal working

CAS #	Analytical Method	Analytical Technique	Sampling Media
8012-95-1	NIOSH 5524	Extraction/GRAV	Pre-weighed PTFE, 37mm
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	768-960	75 µg	
Interferences	Comments		
The method is non-specific and measures all substances extractable by organic solvents.	NIOSH 5524 recommends submitting one bulk sample of each type of fluid for solubility testing. After sampling, ship and store cold. Samples must be analyzed within two weeks of collection.		

Mineral Spirits (Stoddard Solvent)

CAS #	Analytical Method	Analytical Technique	Sampling Media
8052-41-3	NIOSH 1550	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1.3-20	1.4 µg	CS ₂
Interferences	Comments		
	Please send bulk sample. Ship bulk sample separately from air samples.		

Mineral Spirits (Stoddard Solvent)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.03		15-480		2.1 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples. Ship and store at controlled room temperature.		

Molybdenum as Mo

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-98-7	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		0.10 µg
Interferences		Compatibility Code		
		Metals		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Molybdenum as Mo

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7439-98-7	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		1.0 µg
Interferences		Compatibility Code		
		Metals 2		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Morpholine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-91-8	NIOSH S-150	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		10		2.8 µg
Interferences		Compatibility Code		
		dil acid		
Interferences		Comments		

Naphthalene

CAS #	Analytical Method	Analytical Technique	Sampling Media
91-20-3	NIOSH 1501	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	0.47 µg	CS ₂
Interferences		Comments	
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.	

Naphthalene

CAS #	Analytical Method	Analytical Technique	Sampling Media
91-20-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.10	15-480	0.70 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Naphthalene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media
91-20-3	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	200-1000	0.32 µg	PNAs
Interferences		Comments	
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.	

Naproxen Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media
22204-53-1		LC-MS	1.0 µm PTFE Filter
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	200	0.0010 µg	
Interferences		Comments	
		It is critical that the PTFE filter used for sampling must have a pore size of 1.0µm. Contact the lab to determine the minimum air volume required. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.	

Naproxen Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media
22204-53-1		LC-MS	GFF Wipes
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
N/A	N/A	0.0020 µg	
Interferences	Comments		
	Follow "Naproxen Sodium Wipe Sampling Procedure. Order media one week ahead of survey. Media are prepared when ordered. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Nickel and inorganic compounds as Ni

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-02-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	50-1000	0.10 µg	Metals
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Nickel and inorganic compounds as Ni

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-02-0	OSHA ID-121 OSHA ID-125G	ICP	ghost wipe
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
NA	NA	0.50 µg	Metals2
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Nicotine

CAS #	Analytical Method	Analytical Technique	Sampling Media
54-11-5	NIOSH 2544	HPLC	XAD-2 (SKC 226-30-04)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	60-400	2.0 µg	
Interferences	Comments		

Nitric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7697-37-2	NIOSH 7903	IC	MCE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-0.5	5-100	2.3 µg	Acid1
Interferences		Comments	
Particulate salts of the acid will give a positive interference.		Use the maximum flow rate for STEL sampling.	

Nitric acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7697-37-2	NIOSH 7907	IC	SKC 225-9032
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	35-600	3.8 µg	Acid3
Interferences		Comments	
Inorganic acids can react with co-sampled particulate matter on the pre-filter, leading to low results. Potentially interfering particulate chlorides and nitrates removed by the pre-filter can react with the sampled acids and liberate HCl and HNO ₃ , which gets collected on the sampling filter, leading to high results.		Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.	

Nitric Oxide and Nitrogen Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
10102-43-9	OSHA ID-190 (NO), OSHA ID-182 (NO ₂)	IC	SKC 226-40A or SKC 226-182A
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	3-24	0.78 µg	NO&NO ₂
Interferences		Comments	
		Use SKC 226-40A (2TEA coated tubes + oxidizer) to sample NO only or NO ₂ and NO. Store and ship cold.	

Nitroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-24-3	NIOSH 2526	GC-FID	XAD-2 tube (SKC 226-3002A) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	1.5-3	1.8 µg	Ethyl Acetate
Interferences		Comments	
		Sample using 2 XAD-2, (front 600 mg and backup 300 mg), tubes in series. After sampling, separate and cap sorbent tubes.	

Nitrogen Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
10102-44-0	OSHA ID-182	IC	SKC 226-40-2
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2	35	1.2 µg	
Interferences		Comments	
		Sample with TEA IMS (226-40-02) for NO ₂ only. Store and ship cold.	

Nitromethane

CAS #	Analytical Method	Analytical Technique	Sampling Media
75-52-5	NIOSH 2527	GC-FID	CS 106 (SKC 226-111A)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	5-10	2.3 µg	Ethyl Actate
Interferences		Comments	

Nitrous Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
10024-97-2	OSHA ID-166	GC-ECD	AT575 (Assay Technology N ₂ O Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.75	15-480	0.20 µg	
Interferences		Comments	
Halogenated anesthetic gases, CFC's and HCFC's do not interfere when tested at their PELs.		Monitors must be received by Lab within one week after sampling and stored at controlled room temperature. Do not use after expiration date. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab.	

Oil Mist (Mineral Oil)

CAS #	Analytical Method	Analytical Technique	Sampling Media
8012-95-1	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		This method is not for oil mist containing PNAs; it does not collect vapor.	

Oil Mist (Mineral Oil) used in metal working

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8012-95-1	NIOSH 5524	Extraction/GRAV	Pre-weighed PTFE, 37mm	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		768-960		0.75 µg
Interferences		Comments		
The method is non-specific and measures all substances extractable by organic solvents.		NIOSH 5524 recommends submitting one bulk sample of each type of fluid for solubility testing. After sampling, ship and store cold. Samples must be analyzed within two weeks of collection.		

Organic Solvent Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	LM-GCMS-13			
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
Interferences		Comments		
		See List of Scans for list of individual organic solvents.		

Oxalic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
144-62-7	OSHA ID-PV2115	IC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.0		30-100		3.0 µg
Interferences		Comments		

Ozone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
10028-15-6	OSHA ID-214	IC	Treated GFF (GFF, NaNO ₂)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.25-0.5 (TWA) 1.5 (STEL)		90-480		4.8 µg
Interferences		Comments		
Particulate salts of nitrate and nitric acid give positive interferences for ozone. SO ₂ will cause a negative interference.		Treated filter is stable for 1 month. Order media at least one week ahead of survey. Media are prepared when ordered. Refrigerate media before and after sampling. Ship and store cold. For 25% of TLV, you will need 190L. Sulfur dioxide (SO ₂) has a negative interference on ozone. If SO ₂ is suspected, check for its presence using a SO ₂ detector tube. If present, request oxidizer tubes in addition to the ozone filters. Oxidizer tube precedes the ozone filter in the sampling train. After sampling seal the cassette with the end plugs and wrap each cassette with foil.		

Palladium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-05-3	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		50-1000	0.32 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Paraffin Wax Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8002-74-2	NIOSH 0500	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15		100-7200	50 µg	
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

Particles (insoluble or poorly soluble) Not otherwise specified; inhalable

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	HSE MDHS 14	GRAV	Pre-weighed PVC, IOM	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2		40-960	100 µg	
Interferences		Comments		
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab 1 week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.		

Particles (insoluble or poorly soluble) Not otherwise specified; respirable

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7		100-816	50 µg	
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Particles (insoluble or poorly soluble) Not otherwise specified; total

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.	

Pentane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
109-66-0	NIOSH 1500	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	4	0.48 µg	CS ₂
Interferences		Comments	

Pentane(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
109-66-0	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
10.8	15-180	0.73 µg	CS ₂
Interferences		Comments	
		Ship and store at controlled room temperature.	

Pentanedione(2,3-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
600-14-6	OSHA 1012	GC-ECD	Silica Gel Tube (SKC 226-183) (2 tubes in series)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	10L (TWA); 3L (short term)	0.040 µg	95% EtOH
Interferences		Comments	
		Samples are collected on two specially washed silica gel tubes in series. Samples should be protected from the light during and after sampling. Separate and cap tubes after sampling. Order aluminum foil for wrapping the samples. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.	

Pentanedione(2,4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-54-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.8	15-300	2.7 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Pentanone(2-) (Methyl propyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-87-9	NIOSH 1300	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	0.70 µg	CS ₂
Interferences	Comments		

Pentanone(2-) (Methyl propyl ketone)

CAS #	Analytical Method	Analytical Technique	Sampling Media
107-87-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.89	15-480	1.1 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Peracetic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-21-0	NON 57	HPLC	SKC 225-9030 (Hydrogen Peroxide) + SKC-226-193-UC (Peracetic Acid)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1 LPM	15 Liters Max	0.50 µg/ sample	Sample with hydrogen peroxide
Interferences	Comments		
Must be collected using SKC 225-9030 pre-filter.	Peracetic acid must be collected using a SKC 225-9030 pre-filter. This analysis is sub-contracted to an outside laboratory. This method is not covered under the laboratory's AIHA LAP, LLC scope of accreditation. Turnaround time is 10 business days.		

Perchloroethylene (Tetrachloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-40		1.8 µg
Interferences		Comments		

Perchloroethylene (Tetrachloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
127-18-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.10		15-480		2.7 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Petroleum Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	NIOSH 1550	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.3-20		0.50 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

Petroleum Ether

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.92		15-480		0.75 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples. Ship and store at controlled room temperature.		

Phenanthrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
85-01-8	OSHA 58	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		0.12 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

Phenanthrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
85-01-8	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.12 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Phenol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-95-2	OSHA 32	HPLC	XAD-7 (SKC 226-95)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		5-24		0.20 µg
Interferences		Comments		
		Phenol and cresol		

Phenylcyclohexene(4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4994-16-5	NIOSH 1500	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-60		0.38 µg
Interferences		Comments		
		Please indicate if samples are for "LEED" compliance. Minimum sample volume for "LEED" samples is 60 L. Preferred air sampling method.		

Phenylcyclohexene (4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
4994-16-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.28		15-480		0.57 µg
Interferences		Comments		
		AT566 passive monitors cannot be used for "LEED" compliance. Ship and store at controlled room temperature.		

Phenylene(1,3-) diamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-45-2	OSHA 87	HPLC	GFF, Acid	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		10-1000		0.83 µg
Interferences		Comments		
		Sampling media has a short shelf-life so it is not kept in stock. Contact lab one week before sampling to order media.		

Phosphine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7803-51-2	OSHA 1003	ICP	GFF- PE,HgCl ₂ (SKC 225-9018)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		250		3.6 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Minimum air volume required at ¼ of TLV is 250. Sampling media has a shelf-life of only two weeks so it is not kept in stock. Contact lab one week before sampling to order media. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Phosphoric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7664-38-2	NIOSH 7903	IC	MCE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.2-0.5		25-100		2.3 µg
Interferences		Comments		
Particulate salts of the acid will give a positive interference.		Sample at a flow rate of 0.5 lpm for STEL.		

Phosphoric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-38-2	NIOSH 7908	IC	SKC 225-9033
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-5	50-1000	3.8 µg	Acid4
Interferences		Comments	
Particulate salts of phosphate will give a positive interference.		Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.	

Phosphorus (elements)

CAS #	Analytical Method	Analytical Technique	Sampling Media
7723-14-0	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	330-1000	3.3 µg	Metals
Interferences		Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Phthalic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media
85-44-9	OSHA 90	HPLC	GFF, Vamine
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.0	75	1.0 µg	
Interferences		Comments	
Isocyanates, acid chlorides and other anhydrides will give a positive interference.		Sample open-faced. Order filters one week ahead, filters are prepared when ordered and have one month shelf-life.	

Piperazine

CAS #	Analytical Method	Analytical Technique	Sampling Media
110-85-0	OSHA In-house IMIS P250	HPLC	XAD-2, NITC (SKC 226-30-18)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1	10	0.12 µg	
Interferences		Comments	
		Turnaround time is 10 business days.	

Platinum Metal and Soluble Salts as Pt

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-06-4	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		500-1000		0.25 µg
Interferences			Compatibility Code	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			Metals	
Interferences			Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

PNA Scan (NIOSH 5506)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		
Interferences			Compatibility Code	
Asphalt fumes will interfere.				
Interferences			Comments	
Asphalt fumes will interfere.			After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold. Analysis price is for the filter and tube together.	

PNA Scan (OSHA 58)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	OSHA 58	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		
Interferences			Compatibility Code	
Asphalt fumes will interfere.				
Interferences			Comments	
Asphalt fumes will interfere.			After sampling, cap and wrap in aluminum foil. Ship and store cold.	

Polychlorobiphenyl (Chlorodiphenyl, 54% Chlorine) (PCB)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
11097-69-1	NIOSH 5503	GC-MS	Florisil Tube & GFF (SKC 226-39 with Swinnex Filter)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-50		1.1 µg
Interferences			Compatibility Code	
Other chlorinated pesticides may interfere in the quantification of PCB.				
Interferences			Comments	
Other chlorinated pesticides may interfere in the quantification of PCB.			This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.	

Polychlorobiphenyl (Chlorodiphenyl, 42% Chlorine) (PCB)

CAS #	Analytical Method	Analytical Technique	Sampling Media
53469-21-9	NIOSH 5503	GC-MS	Florisil Tube & GFF (SKC 226-39 with Swinnex Filter)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.2	1-50	0.97 µg	
Interferences	Comments		
Other chlorinated pesticides may interfere in the quantification of PCB.	This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Polymeric (Methylene Bisphenyl Isocyanate) (Polymeric MDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media
9016-87-9	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	15	0.10 µg	Isocyanate
Interferences	Comments		
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.	Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.		

Polyvinyl Chloride (PVC)

CAS #	Analytical Method	Analytical Technique	Sampling Media
9002-86-2	NIOSH 0600	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7	100-816	50 µg	
Interferences	Comments		
All other respirable dusts will interfere.	Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Portland Cement

CAS #	Analytical Method	Analytical Technique	Sampling Media
65997-15-1	NIOSH 0600	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7 (See comment)	500-816	50 µg	
Interferences	Comments		
All other dusts will interfere.	Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Potassium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1310-58-3	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		85-100		17 µg
Interferences		Comments		
All forms of potassium are quantified. Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Propanol(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	NIOSH 1401	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.83 µg
Interferences		Comments		
		Store and ship cold overnight.		

Propanol(n-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.8		15-420		1.2 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Propionaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
123-38-6	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.58		15-480		0.029 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Propionaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-38-6	NIOSH 2016	HPLC	Sep-Pak (WAT047205)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	10-100	0.15 µg	Aldehyde
Interferences		Comments	
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling.	

Propionaldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media
123-38-6	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.1-1.5	1-15	0.059 µg	Aldehyde
Interferences		Comments	
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling.	

Propionic Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-09-4	NIOSH 2011	IC	PTFE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05-0.5	15-100	2.3 µg	Acid2
Interferences		Comments	
		Do not sample with inorganic acids.	

Propoxyethanol(2-) (Ethylene glycol monopropyl ether)

CAS #	Analytical Method	Analytical Technique	Sampling Media
2807-30-9	NIOSH 1403	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.05	1-10	1.1 µg	5%MeOH/MC
Interferences		Comments	

Propyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	OSHA PV2061	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		12		1.3 µg
Interferences		Comments		
		1%DMF/CS ₂		

Propyl Bromide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
106-94-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.26		15-480		2.0 µg
Interferences		Comments		
		CS ₂		
		Ship and store at controlled room temperature.		

Propyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-60-4	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		1-10		0.78 µg
Interferences		Comments		
		CS ₂		
		Preferred for STEL sampling. Use a flow rate of 0.2 lpm for STEL.		

Propyl(n-) Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
109-60-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.08		15-480		1.2 µg
Interferences		Comments		
		CS ₂		
		Ship and store at controlled room temperature.		

Propyl(n-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	NIOSH 1401	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.83 µg
Interferences		Comments		
		1%IPA/CS ₂		
		Store and ship cold overnight.		

Propyl(n-) Alcohol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
71-23-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.8		15-420		1.2 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	NIOSH 1403	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-10		1.3 µg
Interferences		Compatibility Code		
		5%MeOH/MC		
Interferences		Comments		

Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
107-98-2	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.66		15-480		2.0 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Propylene Glycol Monomethyl Ether Acetate (PGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-65-6	NIOSH 1450	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		1.1 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		

Propylene Glycol Monomethyl Ether Acetate (PGMEA)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-65-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.98		15-480		1.7 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Propylene Glycol (1,2-Propanediol)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
57-55-6	NIOSH 5523	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-2		5-60		6.2 µg
Interferences		Comments		

Propylene Oxide (1,2-Epoxypropane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-56-9	ASTM D5578-04	GC-FID	ORBO 78 (SUPELCO 20355)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.15		1-24		0.51 µg
Interferences		Comments		
		Sample separately from CS ₂ compatible solvents. Store and ship cold overnight.		

Propylene Oxide (1,2-Epoxypropane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-56-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
12.1		15-480		0.76 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
129-00-0	OSHA 58	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		960		0.60 µg
Interferences		Comments		
Asphalt fumes will interfere.		After sampling, cap and wrap in aluminum foil. Ship and store cold.		

Pyrene (see PNA scan)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
129-00-0	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200-1000		0.60 µg
Interferences		Compatibility Code		
Asphalt fumes will interfere.		PNAs		
Interferences		Comments		
		After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold.		

Pyrethrum

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8003-34-7	NIOSH 5008	HPLC	Glass Fiber Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		20-400		0.13 µg
Interferences		Compatibility Code		
Interferences		Comments		

Pyridine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-86-1	NIOSH 1613	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-1.0		18-150		0.23 µg
Interferences		Compatibility Code		
		MC		
Interferences		Comments		

Resin Acids

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8050-09-7	LM-LC-30	HPLC	GFF, IOM	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
2		200 (minimum)		0.091 µg
Interferences		Compatibility Code		
Interferences		Comments		
		Resin Acids include Abietic Acid and Dehydroabietic Acid. 200-L will give you 46% (0.00046 mg/cu m) of TLV. 400-L will give you 23% (0.00023 mg/cu m) of TLV. 900-L will give you 10% (0.00010 mg/cu m) of TLV.		

Resorcinol

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-46-3	OSHA PV2053	GC-FID	OVS 7 (SKC 226-57)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1		20-400		5.0 µg
Interferences			Compatibility Code	
			MeOH	
Interferences			Comments	

Rhodium as Rh

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-16-6	OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		500-1000		0.50 µg
Interferences			Compatibility Code	
			Metals	
Interferences			Comments	
Spectral interferences are the primary interferences encountered in ICP-AES analysis.			As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.	

Scan for Aldehydes

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 2016	HPLC		
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
Interferences			Compatibility Code	
Interferences			Comments	
			See List of Scans for individual aldehydes.	

Scan for Inorganic Acids

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 7903	IC	MCE/SKC 226-10-03	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.2-0.5		50-100		
Interferences			Compatibility Code	
Interferences			Comments	
			See List of Scans for individual inorganic acids/anions.	

Scan for Isocyanates

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	15		
Interferences	Comments		
	Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold. See List of Scans at the back of the guide for list of individual isocyanate.		

Scan for Organic Solvents

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	LM-GCMS-13	GC-FID GC-MS	
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
Interferences	Comments		
	See List of Scans for individual organic solvents.		

Scan for PNAs (NIOSH 5506)

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 5506	HPLC	PTFE/XAD-2 Tube (PTFE/SKC 226-30-04)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	200-1000		
Interferences	Comments		
Asphalt fumes will interfere.	After sampling, separate filter from sorbent tube. Cap and wrap individually in aluminum foil. Ship and store cold. Price is for the analysis of filter and tube together.		

Scan for PNAs (OSHA 58)

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	OSHA 58	HPLC	Glass Fiber Filter
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	960		
Interferences	Comments		
Asphalt fumes will interfere.	After sampling, cap and wrap in aluminum foil. Ship and store cold.		

Selenium and Compounds as Se

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7782-49-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		120-1000		1.0 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Sevoflurane (Sevofrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28523-86-6	OSHA 106	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05		3		2.9 µg
Interferences		Comments		
		Store and ship cold overnight.		

Sevoflurane (Sevofrane)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
28523-86-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.89		15-480		4.4 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Silica Cristobalite

CAS #	Analytical Method	Analytical Technique	Sampling Media	
14464-46-1	NIOSH 7500	XRD	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
see comments		600-1200		7.5 µg
Interferences		Comments		
Mica, potash, feldspars, zircon, graphite and aluminosilicates will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm; 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm. Dusty atmospheres require much smaller sample volumes (<600 L) to obtain optimum dust loading of 2 mg on filter. 600-L will give you 52% (0.013 mg/cu m) of TLV and 900-L for 33% (0.0083 mg/cu m) of TLV. Sample 600-L to quantify at 25% of the OSHA PEL.		

Silica Quartz

CAS #	Analytical Method	Analytical Technique	Sampling Media
14808-60-7	NIOSH 7500	XRD	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
see comments	600-1200	7.5 µg	Silica
Interferences	Comments		
Mica, potash, feldspars, zircon, graphite and aluminosilicates will interfere.	Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm; 3-piece cassette for BMRC (SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm. Dusty atmospheres require much smaller sample volumes (<600 L) to obtain optimum dust loading of 2 mg on filter. Sample 600-L for 52% (0.013 mg/cu m) of the TLV and 900-L for 33% (0.0083 mg/cu m) of TLV. Sample 600-L to quantify at 25% of the OSHA PEL.		

Silver Metal and Soluble Compounds as Ag

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-22-4	NIOSH 7300	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4	250-1000	0.50 µg	Metals1
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Soapstone

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	34-7200	10 µg	
Interferences	Comments		
All other dusts will interfere.	For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

Soapstone

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0600	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1.7	67-816	50 µg	
Interferences	Comments		
All other respirable dusts will interfere.	Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (British Medical Research Council-SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Sodium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-23-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		60-1000		7.5 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Sodium Hydroxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1310-73-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		60-1000		13 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of sodium are quantified. Use a flow rate of 2 lpm for STEL. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Starch

CAS #	Analytical Method	Analytical Technique	Sampling Media	
9005-25-8	NIOSH 0500	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-15		20-7200		50 µg
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 LPM, for area sampling up to 15 LPM.		

Stoddard Solvent

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	NIOSH 1550	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.3-20		2.3 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

Stoddard Solvent

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8052-41-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.03		15-480		3.4 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples. Ship and store at controlled room temperature.		

Strontium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-24-6	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		620-1000		0.052 µg
Interferences		Compatibility Code		
		Metals		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Styrene (Vinyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-42-5	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.02-1		1-10		0.43 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%. Store and ship cold overnight. 2019 NIC, TWA = 10ppm, STEL = 20 ppm, OTO, A3 , BEI adopted in 2020.		

Styrene (Vinyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
100-42-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.21		15-480		0.65 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Sulfur Dioxide

CAS #	Analytical Method	Analytical Technique	Sampling Media
7446-09-5	OSHA ID 1011	IC	SKC 226-177
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.05 for TWA, 0.5 for STEL	12 – 25	1.5 µg	
Interferences	Comments		
Particulate salts of sulfate, sulfur trioxide and sulfuric acid may give positive interferences for sulfur dioxide.	Order media at least one week ahead of time. Specialty Filter. Media charge applies.		

Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-93-9	OSHA ID-113	IC	PPI-MCE Filter (Thoracic)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2.0	90-230	4.6 µg	Acid1
Interferences	Comments		
Particulate salts of sulfate will give a positive interference.	The TLV for sulfuric acid is as thoracic particulate mass. Use 37-mm MCE 3 piece cassette for BGI GK2.69 cyclone at 1.6 lpm or request thoracic PPI pre-loaded with 0.8µm MCE filter. Order media at least one week ahead of time. Specialty Filter. Media charge applies.		

Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-93-9	NIOSH 7903	IC	MCE/SKC 226-10-03
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.2-0.5	50-150	4.6 µg	Acid1
Interferences	Comments		
Particulate salts of sulfate will give a positive interference.			

Sulfuric Acid

CAS #	Analytical Method	Analytical Technique	Sampling Media
7664-93-9	NIOSH 7908	IC	SKC 225-9033
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-5	50-500	4.6 µg	Acid4
Interferences	Comments		
Particulate salts of sulfate will give positive interference.	Order media at least one week ahead of time. Refrigerate media before and after sampling. Ship and store cold. Specialty Filter. Media charge applies.		

Synthetic Vitreous Fibers (Fiber Count)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 7400	PCM	25mm MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.5-16		50-720		0.050 fiber/fld
Interferences		Comments		
All other fibers will cause a positive interference. Chain-like particles may appear fibrous and high levels of non-fibrous dust particles may obscure fibers.		Adjust sampling flow rate and time to obtain optimum fiber loading on the filter. Do not overload filter. Sample open faced. When shipping your samples, do not pack with untreated polystyrene as can lead to fiber loss from electrostatic effect. This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 10 business days.		

Talc

CAS #	Analytical Method	Analytical Technique	Sampling Media	
14807-96-6	NIOSH 0600	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1.7		100-816		50 µg
Interferences		Comments		
All other respirable dusts will interfere.		Use pre-weighed PVC 2-piece cassette for MSA (Dorr-Oliver) cyclones at 1.7 lpm and 3-piece cassette for BMRC (British Medical Research Council-SKC) cyclones at 2.5 lpm and for BGI-4L at 2.2 lpm.		

Tantalum and Tantalum Oxide Dust as Ta

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-25-7	NIOSH 0500	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		100-1000		50 µg
Interferences		Comments		
All other dusts will interfere.				

Tellurium and Compounds as Te

CAS #	Analytical Method	Analytical Technique	Sampling Media	
13494-80-9	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		160-1000		0.54 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Excludes hydrogen telluride. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Tetrachloroethylene (Perchloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
127-18-4	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-40	1.8 µg	CS ₂
Interferences	Comments		

Tetrachloroethylene (Perchloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
127-18-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
7.10	15-480	2.7 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Tetrahydrofuran [THF]

CAS #	Analytical Method	Analytical Technique	Sampling Media
109-99-9	NIOSH 1609	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-9	0.83 µg	CS ₂
Interferences	Comments		
	Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm. High humidity may greatly decrease breakthrough volume.		

Tetrahydrofuran

CAS #	Analytical Method	Analytical Technique	Sampling Media
109-99-9	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
10.8	15-480	1.2 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Thallium and Compounds, as Tl

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-28-0	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	PVC, IOM MCE, IOM	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		500-1000		1.0 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. Contact Lab 1 week before intended use. The availability of IOM samplers is limited. rental charge for IOM samplers applies.		

Thiram

CAS #	Analytical Method	Analytical Technique	Sampling Media	
137-26-8	NIOSH 5005	HPLC	PTFE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		20-400		3.6 µg
Interferences		Comments		
		Ship samples to the lab immediately after sampling, via overnight shipping. Samples are stable for 7 days at 25°C.		

Tin and Compounds as Sn

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-31-5	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		0.52 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Includes metal, metal oxide, inorganic tin compounds (except tin hydride) and organic compounds. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Titanium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-32-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		25-1000	0.25 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		See titanium dioxide entry for special instructions on collection of titanium dioxide. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Titanium Dioxide

Titanium Dioxide				
CAS #	Analytical Method	Analytical Technique	Sampling Media	
13463-67-7	OSHA ID-121	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		25-1000	0.63 µg	
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of titanium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples. Sample separately from other metals.		

Toluene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-88-3	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2		1-10	0.44 µg	CS ₂
Interferences		Comments		
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.		

Toluene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-88-3	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
9.48		15-480	0.67 µg	CS ₂
Interferences		Comments		
		Ship and store at controlled room temperature.		

Toluene-2,4-diisocyanate (2,4-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media
584-84-9	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	15	0.010 µg	Isocyanate
Interferences		Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.	

Toluene-2,4- diisocyanate (2,4-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media
584-84-9	OSHA 5002	HPLC	Glass Fiber Filter Wipes
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
N/A	N/A	0.010 µg	Isocyanate
Interferences		Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead of survey. Media are prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.	

Toluene-2,6-diisocyanate (2,6-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media
91-08-7	OSHA 5002	HPLC	Isocyanate Filter (GFF, 1-2PP)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1	15	0.010 µg	Isocyanate
Interferences		Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Keep media refrigerated before and after sampling. Sample open-faced. Ship cold overnight. Filter is stable for 6 months from date prepared if kept cold.	

Toluene-2,6- diisocyanate (2,6-TDI)

CAS #	Analytical Method	Analytical Technique	Sampling Media
91-08-7	OSHA 5002	HPLC	Glass Fiber Filter Wipes
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
N/A	N/A	0.010 µg	Isocyanate
Interferences		Comments	
Potential interferences include anhydrides, amines, alcohols and carboxylic acids.		Follow "Isocyanate Wipe Sampling Procedure". Immediately after sampling, glass fiber filters must be placed in a vial containing derivatizing solution. Order media one week ahead, media is prepared when ordered. Derivatizing solution has a shelf life of 6 months if kept cold.	

Toluidine(o-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-53-4	NIOSH 2002	GC-FID	SGT (SKC 226-10)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.02-0.5		5-30		1.4 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

Tributyl Phosphate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
126-73-8	NIOSH 5034	GC-FID	MCE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-3		2-100		0.41 µg
Interferences		Comments		
		Sample separately from CS ₂ compatible solvents.		

Trichloro(1,1,2-)-1,2,2-trifluoroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
76-13-1	NIOSH 1020	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.05		1-10		2.9 µg
Interferences		Comments		

Trichloro(1,1,2-)-1,2,2-trifluoroethane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
76-13-1	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
6.69		15-480		4.4 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Trichlorobenzene(1,2,4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
120-82-1	NIOSH 1003	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1-10		0.29 µg
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Trichloroethane(1,1,1-) (Methyl Chloroform)

CAS #	Analytical Method	Analytical Technique	Sampling Media
71-55-6	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	2.2 µg	CS ₂
Interferences	Comments		
	Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Trichloroethane(1,1,1-) (Methyl Chloroform)

CAS #	Analytical Method	Analytical Technique	Sampling Media
71-55-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
7.93	15-480	3.3 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Trichloroethane(1,1,2-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-00-5	NIOSH 1003	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	2-60	1.6 µg	CS ₂
Interferences	Comments		

Trichloroethane(1,1,2-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-00-5	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
7.58	15-480	2.9 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Trichloroethylene

CAS #	Analytical Method	Analytical Technique	Sampling Media
79-01-6	NIOSH 1022	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	2.0 µg	CS ₂
Interferences	Comments		

Trichloroethylene

CAS #	Analytical Method	Analytical Technique	Sampling Media	
79-01-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.60		15-480		3.0 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Triethanolamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
102-71-6	NIOSH 2007	IC	Silica Gel Tube (SKC 226-10-03)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.5		15-25		7.5 µg
Interferences		Comments		
		Store in freezer after sampling. Ship and store cold.		

Triethylamine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
121-44-8	OSHA PV2060	GC-FID	XAD-7, Acid (SKC 226-98)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		10		4.0 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.		This analysis is sub-contracted to an AIHA LAP, LLC accredited lab. Standard turnaround time is 15 business days.		

Triethylenetetramine

CAS #	Analytical Method	Analytical Technique	Sampling Media	
112-24-3	OSHA 60	HPLC	XAD-2, NITC (SKC 226-30-18)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1		10		0.24 µg
Interferences		Comments		
Nitrogen compounds that co-elute will interfere.				

Triglycidyl Isocyanurate(1,3,5)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
2451-62-9	Ciba-Geigy C321A	GC-MS	PTFE Filter	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1		100-240	4.0 µg	Acetone
Interferences		Comments		
		Samples can be collected using either 25mm or 37mm PTFE filters. Turnaround is 10 business days.		

Trimellitic Anhydride

CAS #	Analytical Method	Analytical Technique	Sampling Media	
552-30-7	OSHA 98	HPLC	GFF, Vamine	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2		480	0.50 µg	
Interferences		Comments		
		Order media one week ahead of survey. Media are prepared when ordered. Sample open-faced.		

Trimethylbenzene(1,2,4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-63-6	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2		1-10	0.44 µg	CS ₂
Interferences		Comments		

Trimethylbenzene(1,2,4-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
95-63-6	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.36		15-480	0.66 µg	CS ₂
Interferences		Comments		
		Ship and store at controlled room temperature.		

Trimethylbenzene(1,3,5-)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-67-8	NIOSH 1501	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2		1-10	0.45 µg	CS ₂
Interferences		Comments		

Trimethylbenzene(1,3,5-)

CAS #	Analytical Method	Analytical Technique	Sampling Media
108-67-8	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
8.36	15-480	0.68 µg	CS ₂
Interferences	Comments		
	Ship and store at controlled room temperature.		

Tungsten and Compounds as W (in the absence of Cobalt)

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-33-7	OSHA ID-213	ICP	MCE Filter
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	30-480	0.25 µg	
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	Recommended air volume for STEL is 30L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Tungsten, as W Soluble Compounds

CAS #	Analytical Method	Analytical Technique	Sampling Media
7440-33-7	OSHA ID-213	ICP	MCE or PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	30-480	0.25 µg	
Interferences	Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.	Recommended air volume for STEL is 30L. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

TVOC as n-Hexane

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 1500	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-10	0.40 µg	CS ₂
Interferences	Comments		
	For IAQ and LEED sampling where TVOC and formaldehyde are collected, do not sample with Charcoal Tubes and DNPH tubes in tandem. The sorbent in the DNPH tubes may off-gas acetonitrile which can cause a positive interference in the TVOC results.		

TVOC as n-Hexane

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.88		15-480		0.60 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

Valeraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	AT571 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.21		15-480		0.043 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Ship and store at controlled room temperature.		

Valeraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	Sep-Pak (WAT047205)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		10-100		0.22 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling. Sample at 1.5 lpm for STEL.		

Valeraldehyde

CAS #	Analytical Method	Analytical Technique	Sampling Media	
110-62-3	NIOSH 2016	HPLC	SGT, DNPH (SKC 226-119)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-1.5		1-15		0.087 µg
Interferences		Comments		
Other aldehydes and ketones will react with the 2,4-DNPH but can be chromatographically resolved.		Keep media refrigerated before and after sampling. Ship cold overnight. Preferred for STEL sampling. Sample at 1.5 lpm for STEL.		

Vanadium Pentoxide as V

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-62-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		120-1000		0.10 µg
Interferences		Compatibility Code		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Metals		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		All forms of vanadium are quantified. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Vanadium

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-62-2	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		0.22 µg
Interferences		Compatibility Code		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		Metals 2		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Vegetable Oil Mist

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 0500	GRAV	Pre-weighed PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-15		20-7200		50 µg
Interferences		Compatibility Code		
All other dusts will interfere.		50 µg		
Interferences		Comments		
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.		

Vinyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-05-4	NIOSH 1453	GC-FID	ORBO 92 Tube	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.1-0.2		0.75-24		1.2 µg
Interferences		Compatibility Code		
		5%MeOH/MC		
Interferences		Comments		
		Preferred for STEL sampling. Sample at a flow rate of 0.2 lpm.		

Vinyl Acetate

CAS #	Analytical Method	Analytical Technique	Sampling Media	
108-05-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.89		15-480		1.8 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Ship and store at controlled room temperature.		

Vinyl Chloride (Chloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-01-4	NIOSH 1007	GC-FID GC-MS	Charcoal Tube (small) (2 Tubes in series)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05		0.7-5		0.29 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comments		
		Sample using 2 charcoal tubes in series. Separate and cap tubes before shipping. Possible loss of sample after 2 or more weeks of storage at room temperature. Ship and store cold immediately.		

Vinyl Chloride (Chloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-01-4	LM-GC-68	GC-FID GC-MS	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
11.6		15-480		0.43 µg
Interferences		Compatibility Code		
		CS ₂		
Interferences		Comment		
		Ship and store at controlled room temperature.		

Vinyl(1-)-2-pyrrolidinone

CAS #	Analytical Method	Analytical Technique	Sampling Media	
88-12-0	NIOSH 1302	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.05-0.2		3-125		0.39 µg
Interferences		Compatibility Code		
		5%MeOH/MC		
Interferences		Comments		

Vinylidene Chloride (1,1-Dichloroethylene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
75-35-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
9.32		15-480		0.60 µg
Interferences		Comments		
		Ship and store at controlled room temperature.		

VM & P Naphtha

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	NIOSH 1550	GC-FID	Charcoal Tube (small, large)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
0.01-0.2		1.3-28		1.8 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples.		

VM & P Naphtha

CAS #	Analytical Method	Analytical Technique	Sampling Media	
8032-32-4	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
8.92		15-480		2.7 µg
Interferences		Comments		
		Please send bulk sample. Ship bulk sample separately from air samples. Ship and store at controlled room temperature.		

Welding Fume Scan

CAS #	Analytical Method	Analytical Technique	Sampling Media	
n/a	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		400-1000		
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		See List of Scans for individual metals. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Welding Fumes, Total

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.	

Wood Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	NIOSH 0500	GRAV	Pre-weighed PVC
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-15	40-7200	50 µg	
Interferences		Comments	
All other dusts will interfere.		For personal sampling use a flow rate of 1-2 lpm, for area sampling up to 15 lpm.	

Wood Dust

CAS #	Analytical Method	Analytical Technique	Sampling Media
n/a	HSE MDHS-14	GRAV	PVC, IOM
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
2	960	100 µg	
Interferences		Comments	
All other dusts will interfere.		Use IOM sampler with pre-weighed PVC. Contact Lab one week before intended use. The availability of IOM samplers is limited. Rental charge for the IOM samplers applies.	

Xylene (Dimethyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media
1330-20-7	NIOSH 1501	GC-FID	Charcoal Tube (small, large)
Sampling Rate (LPM) (CC/Min)	Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
0.01-0.2	1-30	1.0 µg	CS ₂
Interferences		Comments	
		Under conditions of high humidity, the breakthrough volumes may be reduced by as much as 50%.	

Xylene (Dimethyl benzene)

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1330-20-7	LM-GC-68	GC-FID	AT566 (Assay Technology Passive Monitor)	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
7.54		15-480		1.5 µg
Interferences		Comments		
		Ship and store at controlled room temperature		

Yttrium and compounds, as Y

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-65-5	NIOSH 7301 NIOSH 7303	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		25-1000		0.025 µg
Interferences		Comments		
		Spectral interferences are the primary interferences encountered in ICP-AES analysis.		

Zinc

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-66-6	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
1-4		120-1000		3.0 µg
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Zinc

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7440-66-6	OSHA ID-121	ICP	Ghost wipe	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)		LOQ
NA		NA		160 ug
Interferences		Comments		
Spectral interferences are the primary interferences encountered in the ICP-AES analysis.		As part of the Lab's QC protocol, yttrium is used as internal standard in ICP analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Zinc Chloride Fume

CAS #	Analytical Method	Analytical Technique	Sampling Media	
7646-85-7	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		120-1000	6.3 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Analysis is for water soluble zinc compounds. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

Zinc Oxide

CAS #	Analytical Method	Analytical Technique	Sampling Media	
1314-13-2	NIOSH 7301 NIOSH 7303 OSHA ID-125G	ICP	MCE or PVC	
Sampling Rate (LPM) (CC/Min)		Sampling Volume (L) (Minutes)	LOQ	Compatibility Code
1-4		120-1000	3.7 µg	Metals
Interferences		Comments		
Spectral interferences are the primary interferences encountered in ICP-AES analysis.		Samples are analyzed for zinc and calculated as zinc oxide. Current TLV is for respirable sampling. As part of the Lab's QC protocol, yttrium is used as internal standard in metal analysis. Please indicate in your sample submission form if yttrium is present in the area where you collected your samples.		

List of analytes by CAS

50-00-0	Formaldehyde	75	60-29-7	Ethyl Ether	71
50-00-0	Formaldehyde	75	62-53-3	Aniline	25
50-00-0	Formaldehyde	75	63-25-2	Carbaryl (SEVIN)	46
50-21-5	Lactic Acid	91	64-17-5	Ethyl Alcohol (Ethanol)	70
50-32-8	Benzo[a]pyrene (see PNA scan)	30	64-17-5	Ethyl Alcohol (Ethanol)	70
50-32-8	Benzo[a]pyrene (see PNA scan)	30	64-18-6	Formic Acid	76
53-70-3	Dibenzo[a,h]anthracene (see PNA scan)	58	64-19-7	Acetic Acid	19
54-11-5	Nicotine	112	64-19-7	Acetic Acid	19
56-23-5	Carbon Tetrachloride (tetrachloromethane)	47	64-67-5	Diethyl Sulfate	63
56-23-5	Carbon Tetrachloride (tetrachloromethane)	47	67-56-1	Methanol (Methyl alcohol)	97
56-55-3	Benz[a]anthracene (see PNA scan)	28	67-56-1	Methanol (Methyl Alcohol)	97
57-55-6	Propylene Glycol (1,2-Propanediol)	129	67-56-1	Methanol (Methyl Alcohol)	97
60-29-7	Ethyl Ether	71	67-56-1	Methyl Alcohol (Methanol)	100
			67-56-1	Methyl Alcohol (Methanol)	100
			67-56-1	Methyl Alcohol (Methanol)	101

67-63-0	Isopropyl Alcohol (Isopropanol)	90	75-18-3	Dimethyl Sulfide	64
67-63-0	Isopropyl Alcohol (Isopropanol)	90	75-2-8	Ethylene Oxide	73
67-64-1	Acetone	19	75-21-8	Ethylene Oxide	73
67-64-1	Acetone	20	75-25-2	Bromoform	34
67-66-3	Chloroform (Trichloromethane)	50	75-34-3	Dichloroethane(1,1)	60
67-66-3	Chloroform (Trichloromethane)	50	75-34-3	Dichloroethane(1,1)	60
68-12-2	Dimethylformamide	65	75-35-4	Vinylidene Chloride (1,1-Dichloroethylene)	151
68-12-2	Dimethylformamide	65	75-52-5	Nitromethane	114
71-23-8	Propanol(n-)	125	75-56-9	Propylene Oxide (1,2-Epoxypropane)	129
71-23-8	Propanol(n-)	125	75-56-9	Propylene Oxide (1,2-Epoxypropane)	129
71-23-8	Propyl(n-) Alcohol	127	75-65-0	Butyl(tert-) Alcohol	41
71-23-8	Propyl(n-) Alcohol	128	75-65-0	Butyl(tert-) Alcohol	41
71-36-3	Butyl(n-) Alcohol	39	76-13-1	Trichloro(1,1,2-)-1,2,2- trifluoroethane	143
71-36-3	Butyl(n-) Alcohol	39	76-13-1	Trichloro(1,1,2-)-1,2,2- trifluoroethane	143
71-43-2	Benzene	29	76-22-2	Camphor	45
71-43-2	Benzene	29	76-22-2	Camphor	45
71-55-6	Methyl Chloroform (1,1,1-Trichloroethane)	101	78-59-1	Isophorone	88
71-55-6	Methyl Chloroform (1,1,1-Trichloroethane)	102	78-59-1	Isophorone	89
71-55-6	Trichloroethane(1,1,1-) (Methyl Chloroform)	144	78-83-1	Isobutyl Alcohol	87
71-55-6	Trichloroethane(1,1,1-) (Methyl Chloroform)	144	78-83-1	Isobutyl Alcohol	87
74-90-8	Hydrogen Cyanide	83	78-92-2	Butyl(sec-) Alcohol	40
75-01-4	Vinyl Chloride (Chloroethylene)	150	78-92-2	Butyl(sec-) Alcohol	41
75-01-4	Vinyl Chloride (Chloroethylene)	150	78-93-3	Butanone(2-); (Methyl Ethyl Ketone)	36
75-04-7	Ethylamine	71	78-93-3	Butanone(2-); (Methyl Ethyl Ketone)	36
75-05-8	Acetonitrile	20	78-93-3	Methyl Ethyl Ketone (2-Butanone, MEK)	102
75-05-8	Acetonitrile	20	78-93-3	Methyl Ethyl Ketone (2-Butanone, MEK)	102
75-07-0	Acetaldehyde	18	78-94-4	Methyl Vinyl Ketone	105
75-07-0	Acetaldehyde	18	79-00-5	Trichloroethane(1,1,2-)	144
75-07-0	Acetaldehyde	19	79-00-5	Trichloroethane(1,1,2-)	144
75-09-2	Methylene Chloride (Dichloromethane)	108	79-01-6	Trichloroethylene	144
75-09-2	Methylene Chloride (Dichloromethane)	108	79-01-6	Trichloroethylene	145
75-09-2	Dichloromethane (Methylene chloride)	61	79-06-1	Acrylamide	21
75-09-2	Dichloromethane (Methylene chloride)	61	79-09-4	Propionic Acid	126
75-12-7	Formamide	76	79-10-7	Acrylic Acid	21
75-15-0	Carbon Disulfide	46	79-20-9	Methyl Acetate	99
75-15-0	Carbon Disulfide	47	79-20-9	Methyl Acetate	99

79-21-0	Peracetic Acid	118	100-42-5	Styrene (Vinyl benzene)	136
79-24-3	Nitroethane	113	100-42-5	Styrene (Vinyl benzene)	136
80-05-7	Bisphenol A	33	100-44-7	Benzyl Chloride	31
80-62-6	Methyl Methacrylate	104	100-44-7	Benzyl Chloride	32
83-32-9	Acenaphthene (see PNA scan)	18	100-51-6	Benzyl Alcohol	31
84-66-2	Diethyl Phthalate	63	100-51-6	Benzyl Alcohol	31
84-74-2	Dibutyl Phthalate	59	100-52-7	Benzaldehyde	28
85-01-8	Phenanthrene (see PNA scan)	120	100-52-7	Benzaldehyde	29
85-01-8	Phenanthrene (see PNA scan)	120	100-52-7	Benzaldehyde	29
85-44-9	Phthalic Anhydride	122	100-61-8	Methyl Aniline	101
86-73-7	Fluorene (see PNA scan)	74	101-14-4	Methylene(4,4')-bis(2 chloroaniline) (MOCA)	108
88-12-0	Vinyl(1-)-2-pyrrolidinone	150	101-68-8	Methylene Bisphenyl Isocyanate (MDI)	107
91-08-7	Toluene-2,6-diisocyanate (2,6-TDI)	142	101-68-8	Methylene Bisphenyl Isocyanate (MDI)	107
91-08-7	Toluene-2,6-diisocyanate (2,6-TDI)	142	101-77-9	Methylene(4,4')-Dianiline (MDA)	108
91-20-3	Naphthalene	111	102-71-6	Triethanolamine	145
91-20-3	Naphthalene	111	105-46-4	Butyl(sec-) Acetate	40
91-20-3	Naphthalene (see PNA scan)	111	105-46-4	Butyl(sec-) Acetate	40
92-52-4	Biphenyl (Diphenyl)	32	105-60-2	Caprolactam	45
92-52-4	Diphenyl (Biphenyl)	66	106-46-7	Dichlorobenzene(p-)	59
95-49-8	Chlorotoluene(o-)	51	106-46-7	Dichlorobenzene(p-)	59
95-50-1	Dichlorobenzene(o-)	59	106-48-9	Chlorophenol(p-)	50
95-50-1	Dichlorobenzene(o-)	59	106-89-8	Epichlorohydrin (1-Chloro- 2,3-epoxy propane)	67
95-53-4	Toluidine(o-)	143	106-89-8	Epichlorohydrin (1-Chloro- 2,3-epoxy propane)	67
95-63-6	Trimethylbenzene(1,2,4-)	146	106-94-5	Bromopropane(1-)	35
95-63-6	Trimethylbenzene(1,2,4-)	146	106-94-5	Bromopropane(1-)	35
96-22-0	Diethyl Ketone (3- Pentanone)	62	106-94-5	Propyl Bromide	127
96-22-0	Diethyl Ketone	62	106-94-5	Propyl Bromide	127
96-33-3	Methyl Acrylate	100	106-95-6	Allyl Bromide	22
96-33-3	Methyl Acrylate	100	106-99-0	Butadiene(1,3-)	35
96-37-7	Methyl Cyclopentane	102	106-99-0	Butadiene(1,3-)	35
97-63-2	Ethyl Methacrylate	71	107-02-8	Acrolein	20
98-00-0	Furfuryl Alcohol	76	107-02-8	Acrolein	21
98-01-1	Furfural	76	107-04-0	Bromo(1-)-2-Chloroethane	34
98-82-8	Cumene	55	107-05-1	Allyl Chloride	23
98-82-8	Cumene	55	107-05-1	Allyl Chloride	23
98-83-9	Methyl Styrene(a-)	104	107-06-2	Ethylene Dichloride (1,2-Dichloroethane)	72
98-83-9	Methyl Styrene(a-)	105	107-06-2	Ethylene Dichloride (1,2-Dichloroethane)	72
100-41-4	Ethyl Benzene	70			
100-41-4	Ethyl Benzene	70			

107-07-3	Ethylene Chlorohydrin (2-Chloroethanol)	72	108-83-8	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)	64
107-07-3	Ethylene Chlorohydrin (2-Chloroethanol)	72	108-83-8	Dimethyl(2,6-)-4-heptanone (Diisobutyl ketone)	65
107-13-1	Acrylonitrile (Vinyl Cyanide)	21	108-87-2	Methylcyclohexane	106
107-13-1	Acrylonitrile (Vinyl Cyanide)	22	108-87-2	Methylcyclohexane	106
107-15-3	Ethylenediamine	73	108-88-3	Toluene	141
107-18-6	Allyl Alcohol	22	108-88-3	Toluene	141
107-18-6	Allyl Alcohol	22	108-90-7	Chlorobenzene	49
107-21-1	Ethylene Glycol	72	108-90-7	Chlorobenzene	49
107-41-5	Hexylene Glycol (2-Methyl- 2,4-pentanediol)	82	108-93-0	Cyclohexanol	56
107-87-9	Methyl Propyl Ketone (2-Pentanone)	104	108-93-0	Cyclohexanol	56
107-87-9	Methyl Propyl Ketone (2-Pentanone)	104	108-94-1	Cyclohexanone	56
107-87-9	Pentanone(2-) (Methyl propyl ketone)	118	108-94-1	Cyclohexanone	57
107-87-9	Pentanone(2-) (Methyl propyl ketone)	118	108-95-2	Phenol	120
107-92-6	Butyric Acid	42	109-60-4	Propyl(n-) Acetate	127
107-98-2	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	128	109-60-4	Propyl(n-) Acetate	127
107-98-2	Propylene Glycol Monomethyl Ether (PGME, 1-Methoxy-2-propanol)	128	109-66-0	Pentane(n-)	117
107-98-2	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	98	109-66-0	Pentane(n-)	117
107-98-2	Methoxy(1-)-2-propanol (Propylene glycol monomethyl ether, PGME)	98	109-86-4	Methoxyethanol(2-) (Methyl cellosolve, EGME)	98
108-05-4	Vinyl Acetate	149	109-86-4	Methoxyethanol(2-) (Methyl cellosolve, EGME)	98
108-05-4	Vinyl Acetate	150	109-89-7	Diethylamine	63
108-10-1	Methyl Isobutyl Ketone (MIBK)	103	109-99-9	Tetrahydrofuran [THF]	139
108-10-1	Methyl Isobutyl Ketone (MIBK)	103	109-99-9	Tetrahydrofuran	139
108-21-4	Isopropyl Acetate	89	110-12-3	Methyl Isoamyl Ketone	103
108-21-4	Isopropyl Acetate	90	110-12-3	Methyl Isoamyl Ketone	103
108-31-6	Maleic Anhydride	94	110-19-0	Isobutyl Acetate	86
108-45-2	Phenylene(1,3-) diamine	121	110-19-0	Isobutyl Acetate	86
108-46-3	Resorcinol	131	110-43-0	Heptanone(2-) (Methyl Amyl Ketone)	80
108-65-6	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	128	110-43-0	Heptanone(2-) (Methyl Amyl Ketone)	80
108-65-6	Propylene Glycol Monomethyl Ether Acetate (PGMEA)	128	110-43-0	Methyl Amyl Ketone (2-Heptanone)	101
108-67-8	Trimethylbenzene(1,3,5-)	146	110-43-0	Methyl Amyl Ketone (2-Heptanone)	101
108-67-8	Trimethylbenzene(1,3,5-)	147	110-49-6	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	99
			110-49-6	Methoxyethyl(2-) Acetate (Methyl cellosolve acetate, EGMEA)	99
			110-54-3	Hexane(n-)	81
			110-54-3	Hexane(n-)	82

110-62-3	Valeraldehyde	148	120-82-1	Trichlorobenzene(1,2,4-)	143
110-62-3	Valeraldehyde	148	121-44-8	Triethylamine	145
110-62-3	Valeraldehyde	148	123-31-9	Hydroquinone (Dihydroxybenzene)	84
110-80-5	Cellosolve (2-Ethoxyethanol)	47	123-31-9	Hydroquinone (Dihydroxybenzene)	85
110-80-5	Cellosolve (2-Ethoxyethanol)	48	123-38-6	Propionaldehyde	125
110-80-5	Ethoxyethanol(2-) (Cellosolve)	68	123-38-6	Propionaldehyde	126
110-80-5	Ethoxyethanol(2-) (Cellosolve)	68	123-38-6	Propionaldehyde	126
110-82-7	Cyclohexane	56	123-42-2	Diacetone Alcohol (4-Hydroxy- 4-methyl-2-pentanone)	57
110-82-7	Cyclohexane	56	123-42-2	Diacetone Alcohol (4-Hydroxy- 4-methyl-2-pentanone)	58
110-85-0	Piperazine	122	123-54-6	Pentanedione(2,4-)	118
110-86-1	Pyridine	130	123-72-8	Butyraldehyde(n-)	42
110-91-8	Morpholine	110	123-72-8	Butyraldehyde(n-)	42
111-15-9	Ethoxyethyl(2-) Acetate (Cellosolve acetate)	68	123-72-8	Butyraldehyde(n-)	42
111-30-8	Glutaraldehyde	77	123-86-4	Butyl(n-) Acetate	39
111-30-8	Glutaraldehyde	78	123-86-4	Butyl(n-) Acetate	39
111-30-8	Glutaraldehyde	78	123-91-1	Dioxane(p-)	65
111-40-0	Diethylenetriamine	63	123-91-1	Dioxane(p-)	66
111-42-2	Diethanolamine	62	124-17-4	Butoxyethoxy(2-(2-)) Ethyl Acetate	37
111-76-2	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	36	126-73-8	Tributyl Phosphate	143
111-76-2	Butoxyethanol(2-); (Butyl Cellosolve); (EGBE)	37	126-98-7	Methylacrylonitrile	106
111-76-2	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)	38	126-99-8	Chloroprene(beta-); (2-Chloro-1,3-butadiene)	50
111-76-2	Butyl Cellosolve; (2-Butoxyethanol); (EGBE)	39	126-99-8	Chloroprene(beta-); (2-Chloro-1,3-butadiene)	50
111-77-3	Methoxyethoxy(2-(2-)) Ethanol (Diethylene glycol methyl ether)	98	127-18-4	Perchloroethylene (Tetrachloroethylene)	119
111-96-6	Methoxyethyl(2-) Ether (Diethylene glycol dimethyl ether)	99	127-18-4	Perchloroethylene (Tetrachloroethylene)	119
112-07-2	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)	37	127-18-4	Tetrachloroethylene (Perchloroethylene)	139
112-07-2	Butoxyethyl(2-) Acetate; (Butyl Cellosolve Acetate); (EGBEA)	37	127-18-4	Tetrachloroethylene (Perchloroethylene)	139
112-07-2	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)	38	127-19-5	Dimethyl Acetamide	64
112-07-2	Butyl Cellosolve Acetate; (2-Butoxyethyl acetate); (EGBEA)	38	127-19-5	Dimethyl Acetamide	64
112-24-3	Triethylenetetramine	145	129-00-0	Pyrene (see PNA scan)	129
112-34-5	Butoxyethoxy(2-(2-)) Ethanol	37	129-00-0	Pyrene (see PNA scan)	130
117-84-1	Diocetyl Phthalate	65	137-26-8	Thiram	140
120-12-7	Anthracene (see PNA scan)	25	138-86-3	Limonene(d-)	93
120-12-7	Anthracene (see PNA scan)	25	138-86-3	Limonene(d-)	93
			140-88-5	Ethyl Acrylate	69
			140-88-5	Ethyl Acrylate	70

141-32-3	Butyl Acrylate	38	584-84-9	Toluene-2,4-diisocyanate (2,4-TDI)	142
141-32-3	Butyl Acrylate	38	584-84-9	Toluene-2,4-diisocyanate (2,4-TDI)	142
141-43-5	Ethanolamine (2-Aminoethanol)	68	600-14-6	Pentanedione(2,3-)	117
141-78-6	Ethyl Acetate	69	624-92-0	Dimethyl Disulfide	64
141-78-6	Ethyl Acetate	69	628-63-7	Amyl Acetate	24
142-82-5	Heptane	79	628-63-7	Amyl Acetate	24
142-82-5	Heptane	80	687-47-8	Ethyl Lactate	71
142-96-1	Dibutyl Ether	58	763-69-9	Ethyl 3-ethoxypropionate	69
144-62-7	Oxalic Acid	115	763-69-9	Ethyl 3-ethoxypropionate	69
151-67-7	Halothane (Fluothane)	79	822-06-0	Hexamethylene Diisocyanate (1,6-) (HDI)	81
151-67-7	Halothane (Fluothane)	79	822-06-0	Hexamethylene Diisocyanate (1,6-) (HDI)	81
191-24-2	Benzo[ghi]perylene (see PNA scan)	30	872-50-4	Methyl(1-)-2-pyrrolidinone	106
205-99-2	Benzo[b]fluoranthene (see PNA scan)	30	872-50-4	Methyl(1-)-2-pyrrolidinone	106
206-44-0	Fluoranthene (see PNA scan)	74	1305-62-0	Calcium Hydroxide	44
207-08-9	Benzo[k]fluoranthene (see PNA scan)	31	1305-78-8	Calcium Oxide	44
208-96-8	Acenaphthylene (see PNA scan)	18	1309-37-1	Iron Oxide	86
218-01-9	Chrysene (see PNA scan)	52	1309-37-1	Iron	86
218-01-9	Chrysene (see PNA scan)	53	1309-48-4	Magnesium Oxide	94
287-92-3	Cyclopentane	57	1310-58-3	Potassium Hydroxide	125
302-01-2	Hydrazine	82	1310-73-2	Sodium Hydroxide	135
431-03-8	Butanedione(2,3-); (Butadione(2,3-), Diacetyl, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)	36	1314-13-2	Zinc Oxide	154
431-03-8	Diacetyl (Biacetyl, 2,3-Butadione, 2,3-Butanedione, Diketobutane, Dimethyl Diketone, Dimethylglyoxal)	58	1317-65-3	Calcium Carbonate	44
540-59-0	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)	60	95-48-7	Cresol, all Isomers	55
540-59-0	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	60	106-44-5		
540-59-0	Dichloroethylene(1,2-)(cis); (Acetylene dichloride)	60	108-39-4		
540-59-0	Dichloroethylene(1,2-)(trans); (Acetylene dichloride)	61	1319-77-3		
540-84-1	Isooctane	88	1321-74-0	Divinyl Benzene	67
540-84-1	Isooctane	88	1330-20-7	Xylene (Dimethyl benzene)	152
540-88-5	Butyl(tert-) Acetate	41	1330-20-7	Xylene (Dimethyl benzene)	153
540-88-5	Butyl(tert-) Acetate	41	1332-21-4	Asbestos (bulk), all forms	27
552-30-7	Trimellitic Anhydride	146	1332-21-4	Asbestos, all forms	27
563-80-4	Methyl Isopropyl Ketone	103	1332-58-7	Kaolin	90
			1333-86-4	Carbon Black	46
			1333-86-4	Carbon Black	46
			1344-28-1	Aluminum Oxide	23
			1344-95-2	Calcium Silicate Synthetic Nonfibrous	45
			1634-04-4	Methyl Tert-butyl Ether (MTBE)	105
			1634-04-4	Methyl Tert-butyl ether (MTBE)	105
			1675-54-3	Diglycidyl Ether of Bisphenol A	63

2426-08-6	Butyl(n-) Glycidyl Ether	40	7440-25-7	Tantalum and Tantalum Oxide Dust as Ta	138
2426-08-6	Butyl(n-) Glycidyl Ether	40	7440-28-0	Thallium and Compounds, as Tl	140
2451-62-9	Triglycidyl Isocyanurate(1,3,5)	146	7440-31-5	Tin and Compounds as Sn	140
2499-95-8	Hexyl Acrylate	82	7440-32-6	Titanium	141
2807-30-9	Propoxyethanol(2-) (Ethylene glycol monopropyl ether)	126	7440-33-7	Tungsten and Compounds as W (in the absence of Cobalt)	147
4098-71-9	Isophorone Diisocyanate (IPDI)	89	7440-33-7	Tungsten, as W Soluble Compounds	147
4098-71-9	Isophorone Diisocyanate (IPDI)	89	7440-36-0	Antimony and compounds as Sb	25
4994-16-5	Phenylcyclohexene(4-)	120	7440-38-2	Arsenic and inorganic compounds, as As	26
4994-16-5	Phenylcyclohexene (4-)	121	7440-38-2	Arsenic and inorganic compounds, as As	26
5124-30-1	Methylene Bis(4-cyclohexylisocyanate) (HMDI)	107	7440-39-3	Barium and soluble compounds as Ba	28
7085-85-0	Ethyl 2-cyanoacrylate	68	7440-41-7	Beryllium and compounds as Be	32
7429-90-5	Aluminum Metal and insoluble compounds	23	7440-41-7	Beryllium and Compounds as Be	32
7439-91-0	Lanthanum	91	7440-43-9	Cadmium and compounds as Cd	43
7439-92-1	Lead and Inorganic Compounds as Pb	92	7440-43-9	Cadmium and compounds as Cd	43
7439-92-1	Lead and Inorganic Compounds as Pb	92	7440-47-3	Chromium and Inorganic Compounds as Cr	51
7439-92-1	Lead and Inorganic Compounds as Pb	92	7440-47-3	Chromium and Inorganic compounds as Cr	51
7439-93-1	Lithium Salts	93	7440-47-3	Chromium, Hexavalent compounds as Cr	52
7439-95-4	Magnesium	94	7440-47-3	Chromium, Hexavalent Compounds as Cr	52
7439-96-5	Manganese, Elemental and Inorganic compounds as Mn	95	7440-48-4	Cobalt and Inorganic compounds as Co	54
7439-96-5	Manganese, elemental and Inorganic compounds as Mn	95	7440-48-4	Cobalt and Inorganic compounds as Co	54
7439-97-6	Mercury as Hg (Elemental and inorganic forms)	96	7440-50-8	Copper (Fume, Dusts and Mists) as Cu	54
7439-97-6	Mercury as Hg (Elemental and inorganic forms)	96	7440-50-8	Copper (Fume, Dusts and Mists) as Cu	54
7439-97-6	Mercury as Hg Particulate	96	7440-56-4	Germanium	77
7439-98-7	Molybdenum as Mo	110	7440-57-5	Gold	78
7439-98-7	Molybdenum as Mo	110	7440-62-2	Vanadium Pentoxide as V	149
7440-02-0	Nickel and inorganic compounds as Ni	112	7440-62-2	Vanadium	149
7440-02-0	Nickel and inorganic compounds as Ni	112	7440-65-5	Yttrium and compounds, as Y	153
7440-05-3	Palladium	116	7440-66-6	Zinc	153
7440-06-4	Platinum Metal and Soluble Salts as Pt	123	7440-66-6	Zinc	153
7440-16-6	Rhodium as Rh	131	7440-69-9	Bismuth	33
7440-22-4	Silver Metal and Soluble Compounds as Ag	134	7440-70-2	Calcium	43
7440-23-5	Sodium	135	7440-74-6	Indium and Compounds as In	85
7440-24-6	Strontium	136	7446-09-5	Sulfur Dioxide	137

7646-85-7	Zinc Chloride Fume	154	8052-41-3	Mineral Spirits (Stoddard Solvent)	109
7647-01-0	Hydrogen Chloride	83	8052-41-3	Mineral Spirits (Stoddard Solvent)	110
7647-01-0	Hydrogen Chloride	83	8052-41-3	Stoddard Solvent	135
7553-56-2	Iodine and Iodides as I	85	8052-41-3	Stoddard Solvent	136
7664-38-2	Phosphoric Acid	121	8052-42-4	Asphalt Fume	27
7664-38-2	Phosphoric Acid	122	8052-42-4	Asphalt Fume as Benzene-Soluble Aerosol	28
7664-39-3	Hydrogen Fluoride, as F	84	9002-86-2	Polyvinyl Chloride (PVC)	124
7664-41-7	Ammonia	24	9005-25-8	Starch	135
7664-93-9	Sulfuric Acid	137	9016-87-9	Methylene Bisphenyl Isocyanate, Polymeric (Polymeric MDI)	107
7664-93-9	Sulfuric Acid	137	9016-87-9	Polymeric (Methylene Bisphenyl Isocyanate) (Polymeric MDI)	124
7664-93-9	Sulfuric Acid	137	10024-97-2	Nitrous Oxide	114
7697-37-2	Nitric Acid	113	10028-15-6	Ozone	115
7697-37-2	Nitric acid	113	10035-10-6	Hydrogen Bromide	82
7722-84-1	Hydrogen Peroxide	84	10035-10-6	Hydrogen Bromide	83
7723-14-0	Phosphorus (elements)	122	10049-04-4	Chlorine Dioxide	48
7726-95-6	Bromine	34	10102-43-9	Nitric Oxide and Nitrogen Dioxide	113
7758-97-6	Lead Chromate as Cr(VI)	93	10102-44-0	Nitrogen Dioxide	114
7778-18-9	Calcium Sulfate (Gypsum)	45	11097-69-1	Polychlorobiphenyl (Chlorodiphenyl, 54% Chlorine) (PCB)	123
7782-42-5	Graphite	79	11097-69-1	Chlorodiphenyl (Polychlorobiphenyl, 54% Chlorine)	49
7782-49-2	Selenium and Compounds as Se	133	12001-26-2	Mica	109
7782-50-5	Chlorine	48	12125-02-9	Ammonium Chloride Fume	24
7783-06-4	Hydrogen Sulfide	84	13463-67-7	Titanium Dioxide	141
7784-42-1	Arsine	26	13494-80-9	Tellurium and Compounds as Te	138
7803-51-2	Phosphine	121	13838-16-9	Enflurane (Ethrane)	67
8002-74-2	Paraffin Wax Fume	116	13838-16-9	Enflurane (Ethrane)	67
8003-34-7	Pyrethrum	130	14807-96-6	Talc	138
8006-61-9	Gasoline	77	14464-46-1	Silica Cristobalite	133
8006-61-9	Gasoline	77	14808-60-7	Silica Quartz	134
8008-20-6	Kerosene	91	22204-53-1	Naproxen Sodium	111
8008-20-6	Kerosene	91	22204-53-1	Naproxen Sodium	112
8012-95-1	Mineral Oil (Oil mist)	109	26675-46-7	Forane (Isoflurane)	74
8012-95-1	Mineral Oil, used in metal working	109	26675-46-7	Forane (Isoflurane)	75
8012-95-1	Oil Mist (Mineral Oil)	114	26675-46-7	Isoflurane (Forane)	87
8012-95-1	Oil Mist (Mineral Oil) used in metal working	115	26675-46-7	Isoflurane (Forane)	88
8032-32-4	Petroleum Ether	119			
8032-32-4	Petroleum Ether	119			
8032-32-4	VM & P Naphtha	151			
8032-32-4	VM & P Naphtha	151			
8050-09-7	Resin Acids	130			

28182-81-2	Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)	80	57041-67-5	Desflurane (Suprane)	57
28182-81-2	Hexamethylene Diisocyanate (1,6-) Homopolymer (HDI Homo)	81	57041-67-5	Desflurane (Suprane)	57
28523-86-6	Sevoflurane (Sevofrane)	133	64742-95-6	Aromatic 100	26
28523-86-6	Sevoflurane (Sevofrane)	133	65996-93-2	Coal Tar Pitch Volatiles, as Benzene Soluble Aerosol	53
34590-94-8	Dipropylene Glycol Methyl Ether (DPGME)	66	65997-15-1	Portland Cement	124
34590-94-8	Dipropylene Glycol Methyl Ether (DPGME)	66	88917-22-0	Dipropylene Glycol Methyl Ether Acetate (DPGMEA)	66
53469-21-9	Chlorodiphenyl (Polychlorobiphenyl, 42% Chlorine)	49			
53469-21-9	Polychlorobiphenyl (Chlorodiphenyl, 42% Chlorine) (PCB)	124			



The illustrations, instructions, and principles contained in the material are general in scope and, to the best of our knowledge, current at the time of publication. Our risk control services are advisory only. We assume no responsibility for: managing or controlling customer safety activities, implementing any recommended corrective measures, or identifying all potential hazards.

No attempt has been made to interpret any referenced codes, standards, or regulations. Please refer to the appropriate government authority for interpretation or clarification.

Insurance underwritten by Liberty Mutual Insurance Co. or its affiliates or subsidiaries.

© 2023 Liberty Mutual Insurance, 175 Berkeley Street, Boston, MA 02116.

libertymutualgroup.com/riskcontrolservices

2023 Edition Lab Sampling Guide