

Request for IH sampling media equipment

Industrial Hygiene Laboratory Services
from Liberty Mutual Insurance

Ship attention

E-Mail

Phone number

Order date

Date needed*

Company

Street address

City

State

Zip

Sampling media: Please indicate quantities needed and include blanks in your total quantity.

Qty. Sorbent Tubes

	Charcoal, small
	Charcoal, large
	Silica Gel (SKC 226-10)
	ORBO 34 (hydrogen sulfide)
	ORBO 77 (ammonia)
	ORBO 78 (ETO)
	Sep paks, WAT047205 (aldehydes) ^{1,2}
	SKC 226-10-03 (alkanolamines)
	SKC 226-10-03 & 25mm MCE prefilter (inorganic acids)
	SKC 226-10-03 & PTFE prefilter (organic acids)
	SKC 226-119, DNPH (aldehydes) ^{1,2}
	SKC 226-183 & foil, order by set: 2 tubes in series (diacetyl) ⁴
	SKC 226-193-UC (peracetic acid) ^{1,5} , includes SKC 225-9030 ^{1,4,5} <i>If sampling for Hydrogen Peroxide only, order SKC 225-9030 (under Specialty Filters)</i>
	SKC 226-30-04 with PTFE Filter & foil (PNAs) ⁴
	SKC 226-30-08, order by set: 2 tubes in series (acrylic acid)
	SKC 226-30-18, NITC (EDA, DETA, TETA)
	SKC 226-39 with Swinnex Filter (PCB)
	SKC 226-40-02 (nitrogen dioxide)
	SKC 226-182A (NO/NO2) ¹
	SKC 226-57 (ethylene glycol) ²
	SKC 226-177 (sulfur dioxide) ^{1,2,5}
	SKC 226-82, order by set: 2 tubes in series (methanol)
	SKC 226-95, XAD 7 (phenol)
	SKC 226-98, acid treated XAD 7 (cyanoacrylates)

Other:

--

Wipe sampling kits

Select from pull down menu:

--

(Indicate number of wipes, templates, gloves needed)

--

Qty. Common Filters

Disposable IOM-PVC ⁵ (Hexavalent Chromium, Inhalable Dust <i>only</i>)		
MCE: 25mm (asbestos)		
MCE: 2-pc, 37mm (metals)		
PVC 25mm, 2-pc pre-weighed		
PVC 2-pc, 37mm pre-weighed		
PVC 3-pc, 37mm pre-weighed		
PPI w/PVC ^{3,5}	Pre-weighed?	
IOM-PVC ^{3,5}	Pre-weighed?	

Qty. Specialty Filters

	GFF with foil (coal tar pitch volatiles)						
	GFF,NaNO2 with foil (ozone) ^{2,4}						
	GFF with 1,2 p,p (isocyanates) ^{1,2}						
	<table><tr><td></td><td>Isocyanate Field Extraction Kit?</td><td></td></tr><tr><td></td><td>Glove Size (for Extraction Kit)?</td><td></td></tr></table>		Isocyanate Field Extraction Kit?			Glove Size (for Extraction Kit)?	
	Isocyanate Field Extraction Kit?						
	Glove Size (for Extraction Kit)?						
	IOM-GFF (resin acids) ⁵	Calibrator?					
	PPI-MCE Thoracic (sulfuric acid) ⁵						
	PTFE, 37mm pre-weighed (metal working fluid/asphalt fume)						
	PTFE/Silver membrane (chlorine & bromine)						
	SKC 225-9030 w/foil (hydrogen peroxide <i>only</i>) ^{1,4,5}						
	SKC 225-9031 (fluorides, aerosol & gases) ^{1,2,5}						
	SKC 225-9032 (hydrogen chloride/bromide, nitric) ^{1,2,5}						
	SKC 225-9033 (phosphoric, sulfuric) ^{1,2,5}						
	QFF (diesel particulate matter – above ground operation)						
	SKC 226-317 with cyclone & holder # cyclones/holders needed: _____ (diesel particulate matter – underground mines) ⁵						

Qty. Passive Vapor Monitors

	555AT, Assay Technology badge (ETO) ²
	566AT, Assay Technology badge (solvents) ²
	571AT, Assay Technology badge (aldehydes) ²
	X575AT, Assay Technology badge (nitrous oxide) ²
	SKC 575-007, Passive sampler (Methanol – long term sampling)
	SKC 575-002, Passive sampler (Methanol – short term sampling)

Notes

- Media must be kept cold
 - Limited shelf-life
 - Please indicate if you need the filters pre-weighed
 - Media must be wrapped in aluminum foil after sampling
 - Media charges apply
- * UPS Ground shipping is free.



Liberty Mutual
INSURANCE

Qty.	Miscellaneous
	Cooler w/ ice packs
	Fritted Bubblers
	KI solution for chlorine dioxide

Qty.	
	Glass Vials w/ parafilm (liquid bulks)
	Plastic Vials (solid bulks)

Equipment

Requires completed Equipment terms and conditions, and request form (RC6128) prior to shipment.

Gas testing

Choose type:

☐ Breathing air

☐ Pipeline

Qty.

	Cylinders
	47mm Teflon filters
	Fixtures

Pump loan

(Hoses, clips, and chargers included)

	Pumps	Rate (Lpm):	
	Flow-limiters		
	Cyclones (use 2pc cassette, rate: 1.7 Lpm)		

Comments/other: (Qty, Description, Sampling for...)

To send this form as an email attachment, first save this PDF with the completed entries to your desktop.

Next, click this button to open an email addressed to the lab. Attach the saved PDF and submit the email.

Liberty Mutual Industrial Hygiene Lab
 94 South Street, Hopkinton MA 01748
Telephone: (800)230-6263 **Fax:** (508) 435-3575
E-Mail: LMIHLaboratory@LibertyMutual.com

